REF	CONT	REF	CONT
CHO-10100	R: 2 x 50 mL	CHO-10300M	R: 6 x 50 mL
CHO-10300	R: 6 x 50 mL	CHO-10300P	R: 6 x 50 mL
CHO-10600	R: 6 x 100 mL		
CHO-10180A	R: 4 x 45 mL		


### INTENDED USE

For the *in vitro* quantitative determination of Total Cholesterol in serum

### METHODOLOGY

Alain and Roeschlau began using cholesterol esterase and oxidase, in a single reagent to determine total cholesterol in serum. Trinder's color system of peroxidase / phenol/4-aminoantipyrine has been used successfully for some time now. The system's only drawback was the corrosive properties of phenol. The present method utilizes a phenol substitute that performs like phenol but without being corrosive.

### PRINCIPLE

The intensity of the red color produced is directly proportional to the total cholesterol in the sample when read at 500nm.

### REAGENT COMPOSITION

Active Ingredients	Concentrations
PIPES	90 mM
4-Aminoantipyrine	0.40 mM
Cholesterol Oxidase	350 U/L
Cholesterol Esterase	500 U/L
Horseradish Peroxidase	1000 U/L
Phenol	28 mM

pH 7.0 ± 0.2

### PRECAUTIONS

- This reagent is for *in vitro* diagnostic use only.
- Reagent contains sodium azide. Poison. Do not ingest.

### REAGENT PREPARATION

Reagent is supplied ready to use.

### REAGENT STORAGE

Store the reagent at 2-8°C (refrigerated).

The reagent is stable until the expiration date when stored at 2-8°C.

### REAGENT DETERIORATION

The reagent should not be used if:

- The reagent is turbid.
- The working reagent does not meet stated performance parameters.

### SPECIMEN COLLECTION AND STORAGE

- Nonhemolyzed serum is recommended.
- Cholesterol in serum is reported stable for seven days at room temperature (18-25°C) and six months when frozen and properly protected against evaporation.

### INTERFERENCES

Studies to determine the level of interference for hemoglobin, bilirubin, and lipemia were carried out, the following results were obtained:

#### Hemoglobin:

No significant interference (± 10%) from hemoglobin up to 200 mg/dL.

#### Bilirubin:

No significant interference (± 10%) from bilirubin up to 19.8 mg/dL.

#### Lipemia :

No significant interference (± 10%) from lipemia up to 1093 mg/dl measured as triglycerides.

### ADDITIONAL EQUIPMENT REQUIRED BUT NOT PROVIDED:

- A clinical chemistry analyzer capable maintaining constant temperature (37°C), and measuring absorbance at 500±20 nm.
- Iron-free deionized water and related equipment, e.g.: pipettes
- Analyzer specific consumables, e.g.: sample and read cups
- Control and calibrator materials

### PROCEDURE

Wavelength	: 500±20 nm
Working temperature	: 37°C
Optical path	: 1 cm
Assay type	: Endpoint
Direction	: Increasing

	Blank	Standard	Sample
Reagent	1000 µL	1000 µL	1000 µL
Distilled water	10 µL	--	--
Standard	--	10 µL	--
Sample	--	--	10 µL

Mix and then incubate for 5 minutes at 37°C. Measure the absorbance of sample and standard against the reagent blank.

### PROCEDURE NOTES

- The reagent and sample volumes may be altered proportionally to accommodate various instrument requirements.
- Grossly lipemic serums require a "sample blank". Add 0.01ml (10µl) of sample to 1.0mL saline, mix and read the absorbance against water. Subtract this value from the patient absorbance to obtain the corrected reading.

### CALCULATION

(Abs = Absorbance)

$\frac{\text{Abs. of sample}}{\text{Abs. of standard}} \times \text{conc. of std. (mg/dL)} = \text{Cholesterol (mg/dL)}$

### LIMITATIONS

Samples with values exceeding 750 mg/dL should be diluted 1:1 with saline, re-run and result multiplied by two.

### CALIBRATION

Aqueous standards can be used to calibrate the procedure or an appropriate serum calibrator.

### QUALITY CONTROL

The integrity of the reaction should be monitored by use of a two level control with known cholesterol values.

### EXPECTED VALUES

Recommended Range:

Desirable Cholesterol:	<200 mg/dL
Borderline-High Cholesterol:	200 – 239 mg/dL
High Cholesterol:	>240 mg/dL

### PERFORMANCE

#### LINEARITY:

When run as recommended the assay is linear to 750 mg/dL.

#### METHOD COMPARISON:

Studies performed between this procedure and a similar methodology yielded the following results:

Number of samples pairs:	49
Range of samples:	1-502 (mg/dL)
Correlation Coefficient:	0.9961
Slope:	1.015
Intercept:	-3.7 (mg/dL)

#### Precision:

Within Run	Level 1	Level 2
Mean (mg/dL)	148	224
S.D. (mg/dL)	2.2	2.1
C.V. (%)	1.5	1.0
Total	Level 1	Level 2
Mean (mg/dL)	148	224
S.D. (mg/dL)	3.5	6.3
C.V. (%)	2.3	2.8

#### SENSITIVITY:

A calibration factor of approximately 1702 was obtained, which is equivalent to a sensitivity of 0.0006 Δ Abs per mg/dL.

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Bilimsel Tıbbi Ürünler Paz.San.Tic.Ltd.Şti.