REF	CONT	REF	CONT
LDL-30080	R1: 2 x 30 mL R2: 2 x 10 mL	LDL-30240M	R1: 4 x 45 mL R2: 2 x 24 mL
LDL-30200	R1: 3 x 50 mL R2: 1 x 50 mL	LDL-30240P	R1: 4 x 45 mL R2: 2 x 24 mL
LDL-30400	R1: 4 x 75 mL R2: 2 x 50 mL		
LDL-30128A	R1: 4 x 24 mL R2: 4 x 8 mL		

### INTENDED USE

This reagent is intended for the *in vitro* quantitative determination of low-density lipoprotein cholesterol (LDL-C) in human serum.

### METHODOLOGY

The assay is based on a modified polyvinyl sulfonic acid (PVS) and polyethylene-glycol methyl ether (PEGME) coupled classic precipitation method with the improvements in using optimized quantities of PVS/PEGME and selected detergents. LDL, VLDL, and chylomicron (CM) react with PVS and PEGME and the reaction results in inaccessibility of LDL, VLDL and CM by cholesterol oxidase (CHOD) and cholesterol esterase (CHER), whereas HDL reacts with the enzymes. Addition of R2 containing a specific detergent releases LDL from the PVS/PEGME complex. The released LDL reacts with the enzymes to produce H<sub>2</sub>O<sub>2</sub> which is quantified by the Trinder reaction.

### REAGENT COMPOSITION

**Reagent 1:** MES buffer (pH 6.5), Polyvinyl sulfonic acid, 4-aminoantipyrine, Cholesterol esterase, Cholesterol oxidase, Peroxidase, Polyethylene-glycol-methyl ester, MgCl<sub>2</sub>, Detergent, EDTA

**Reagent 2:** MES buffer (pH 6.5), TODB N,N-Bis(4-sulfobutyl)-3- methylaniline, Detergent, EDTA

### REAGENT PREPARATION

Reagent 1 and Reagent 2 are ready to use.

### Precautions and Warnings:

- For *in vitro* diagnostic use only.
- DO NOT pipette by mouth.
- Avoid contact with skin and eyes. If spilt, thoroughly, wash affected areas with water.
- All specimens used in the test should be considered potentially infectious.
- Do not use the reagent after the expiration date printed on the kit.

### REAGENT STORAGE

- Store the reagent at 2-8°C.
- The reagent is stable until the expiration date printed on the label when stored at 2-8°C protected from light and contamination prevented during their use.

### REAGENT DETERIORATION:

Do Not Use Reagent If:

- The reagent develops turbidity.
- The reagent fails to recover control values within the assigned range.

### SPECIMEN COLLECTION AND HANDLING

Use fresh patient serum or plasma samples (EDTA, Citrate). Fasting and non-fasting samples can be used.

Storage: If not analyzed promptly, specimens may be stored at 2-8°C for up to 7 days.

### INTERFERENCES

The following substances normally present in serum produced less than 10% deviation at the listed concentrations: Triglyceride at 1000 mg/dL, Ascorbic Acid at 10 mM, Bilirubin at 40 mg/dL, Bilirubin Conjugate at 30 mg/dL, Hemoglobin at 1000 mg/dL.

### Materials Required but not Provided

- A clinical chemistry analyzer or photometre.
- Analyzer specific consumables, eg: sample cups.
- Control and calibrator materials.

### PROCEDURE

Primary Wavelength	: 600 nm
Secondary Wavelength	: 700 nm
Working temperature	: 37°C
Optical path	: 1 cm
Assay type	: Endpoint/sample blank
Direction	: Increasing

- Bring reagent to room temperature (15-30°C).

	Standard	Sample
Reagent 1	300 µL	300 µL
Standard	4 µL	--
Sample	--	4 µL

- Mix and then incubate for 5 min at 37°C, read the absorbance of the sample (Abs.1 sample) and standard (Abs.1 standard) against distilled water.

	Standard	Sample
Reagent 2	100 µL	100 µL

- Return tube to 37° C and after 5 min. read the final absorbance of the sample (Abs.2 sample ) and standard (Abs.2 standard) against distilled water.

- Determine the Δ Abs. for Sample and Standard.

### CALCULATIONS

(ΔAbs. = Abs.2 – Abs.1)

$\frac{\Delta \text{Abs. Sample}}{\Delta \text{Abs. Standard}} \times \text{Concentration of standard} = \text{LDL-Cholesterol (mg/dL)}$

(mg/dL)

### Limitations

- Anticoagulants containing heparin should not be used.
- Protect the reagents from direct sunlight.
- Store the reagents as per instructions..

### CALIBRATION

**bt products** Chem-Calibrator or HDL/LDL Calibrator is required for calibration. If control results are found to be out of range, the procedure should be recalibrated.

### QUALITY CONTROL

Reliability of test results should be routinely monitored with control materials that reasonably emulate performance of patient specimens. Quality control materials are intended for use only as monitors of accuracy and precision. Controls should be run with every working shift in which LDL-Cholesterol assays are performed.

### RESULTS

LDL-Cholesterol concentration is expressed as mg/dL.

To convert from conventional units to S.I.units, multiply the conventional units by 0.02586.

mg/dL x 0.02586 = mmol/L LDL-Cholesterol

mmol/L x 38.66 = mg/dL

### EXPECTED VALUES

- <100 mg/dL Optimal
- 100 - 129 mg/dL Near optimal / above optimal
- 130 - 159 mg/dL Borderline high
- 160 - 189 mg/dL High
- ≥190 mg/dL Verv high

Each laboratory must establish its own range of expected values.

### PERFORMANCE

#### 1. Linearity: 4 – 250 mg/dL

If samples contain LDL cholesterol greater than 250 mg/dL, they should be diluted with saline and retested.

#### 2. Accuracy:

The performance of this assay was compared with the performance of a legally marketed LDL-Cholesterol assay using serum samples. Seventy nine serum samples ranging from 4.0 - 232.0 mg/dL gave a correlation coefficient of 0.9804. Linear regression analysis gave the following equation:  $y = 1.0883x + 0.6078$

#### 3. Precision:

The precision of the **bt products** LDL-Cholesterol Reagent was evaluated according to Clinical Laboratory Standards Institute EP5-A guideline. In the study, three serum specimens containing 95 mg/dL, 146 mg/dL and 210 mg/dL LDL-Cholesterol were tested on Hitachi 917 with 2 runs per day with duplicates over 20 working days. This method has not been tested or certified by the Cholesterol Reference Method Laboratory Network.

	Level 1 95 mg/dL	Level 1 146 mg/dL	Level 1 210 mg/dL
Within- Run Precision :	LDL	LDL	LDL
N	80	80	80
Mean	97,14	147,37	211,47
SD	1,00	1,19	1,38
CV %	1,0%	0,8%	0,7%

	Level 1 95mg/dL	Level 1 146 mg/dL	Level 1 210 mg/dL
Within- Laboratory Precision :	LDL	LDL	LDL
N	80	80	80
Mean	97,14	147,37	211,47
SD	1,55	2,23	2,98
CV %	1,6%	1,5%	1,4%

### REFERENCES

- "Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult treatment Panel III)", JAMA, 285:2486 (2001).
- Crouse JR et al., Studies of low density lipoprotein molecular weight in human beings with coronary artery disease, J. Lipid Res., 26; 566 (1985).
- Badimon 11, Badimon L., Fuester V., Regression of Atherosclerotic Lesions by High-Density Lipoprotein Plasma Fraction in the Cholesterol-Fed Rabbit, Journal of Clinical Investigation, 85: i 234-4 i (I 990).
- Castelli WP et al., Cholesterol and other lipids in coronary heart disease, Circulation, 55; 767 (1977).
- Barr DP, Rliss EM, Eder HA, Protein-lipid relationships in human plasma, Am. J. Med., 1 1;480 (1951).
- Gordon T. et al., High density lipoprotein as a protective factor against coronary heart disease, Am.J.Med.,62;707 (1977).
- William P., Robinson D., Baily A., High density lipoprotein and coronary risk factor, Lancet, 1;72 (I 979).
- Kannel WB, Castelli WP, Gordon T., Cholesterol in the prediction of atherosclerotic disease; New perspectives based on the Framingham study, Am. Intern. Med.90;85 (1979).
- Hongbing Xiao Method and composition for determining low density lipoprotein cholesterol, Chinese Patent CN1379234A (2002).
- Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of the High Blood Cholesterol in Adults (Adult Treatment Panel III).
- Pisani T, Gebksi CP, Leary Et, et al., Accurate Direct Determination of Low-Density Lipoprotein Cholesterol Assay. Arch Pathol Lab Med 1995; 119:1127.

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### BİLİMSSEL KAYNAKÇA

1. "Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult treatment Panel III)", JAMA, 285:2486 (2001).
2. Crouse JR et al., Studies of low density lipoprotein molecular weight in human beings with coronary artery disease, J. Lipid Res., 26; 566 (1985).
3. Badimon 11, Badimon L., Fuester V., Regression of Atherosclerotic Lesions by High-Density Lipoprotein Plasma Fraction in the Cholesterol-Fed Rabbit, Journal of Clinical Investigation, 85: i 234-4 i (I 990).
4. Castelli W.P et al., Cholesterol and other lipids in coronary heart disease, Circulation, 55; 767 (1977).
5. Barr DP, Rliss EM, Eder HA, Protein-lipid relationships in human plasma, Am. J. Med., 1 1;480 (1951).
6. Gordon T. et al., High density lipoprotein as a protective factor against coronary heart disease, Am.I.Med.,62:707 (1977).
7. William P., Robinson D., Baily A., High density lipoprotein and coronary risk factor, Lancet, 1;72 (I 979).
8. Kannel WB, Castelli WP, Gordon T., Cholesterol in the prediction of atherosclerotic disease; New perspectives based on the Framingham study, Am. Intern. Med.90:85 (1979).
9. Hongbing Xiao Method and composition for determining low density lipoprotein cholesterol, Chinese Patent CN1379234A (2002).
10. Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of the High Blood Cholesterol in Adults (Adult Treatment Panel III).
11. Pisani T, GebSKI CP, Leary Et, et al. Accurate Direct Determination of Low-Density Lipoprotein Cholesterol Assay. Arch Pathol Lab Med 1995; 119:1127.

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Manufacturer / Üretici Bilimsel Tıbbi Ürünler Paz.San.Tic.Ltd.Şti.