REF	CONT	REF	CONT	REF	CONT
ALB-10100	2 x 50 mL	ALB-10180M	4 x 45 mL	ALB-10100A2	5 x 20 mL
ALB-10300	6 x 50 mL	ALB-10180P	4 x 45 mL		
ALB-10600	6 x 100 mL	ALB-10240M2	6 x 40 mL		
ALB-10180A	4 x 45 mL	ALB-10240M3	6 x 40 mL		

KULLANIM AMACI

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BİLİMSEL KAYNAKÇA

1. Tietz, N., Textbook of Clinical Chemistry, Philadelphia, W.B. Saunders, 1986 pp. 701-704
2. Webster D: 177. The Immediate Reaction between Bromocresol Green and Serum as a Measure of Albumin Content. Clin Chem 23:663
3. Dumas BT, Warson WA, and Biggs NG: 1971. Albumin Standards and the measurement of Serum Albumin with Bromocresol Green Clin Chem Acta 31:87
4. Young DS, Effects of Drugs on Clinical Laboratory Tests. Third Edition 1990:12-6

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INTENDED USE

For the in vitro quantitative determination of albumin concentration in serum.

METHODOLOGY

At a controlled pH, bromocresol green forms a colored complex with albumin. The intensity of color at 630 nm is directly proportional to albumin content. The instantaneous initial absorbance is obtained as suggested by Webster

Principle: Controlled pH

BCG+ ALBUMIN $\xrightarrow{\text{Controlled pH}}$ Green BCG/ALBUMIN Complex

REAGENT COMPOSITION

Active Ingredients

	Concentrations
Bromocresol green	0,25 mM
Succinat Buffer	85 mM
Surfactant	
pH 4.20 ± 0.1	

PRECAUTIONS

- For *in vitro* diagnostic use only.
- DO NOT pipette by mouth. Avoid contact with skin and eyes. If spilt, thoroughly, wash affected areas with water.
- Reagent contains Sodium Azide as a preservative. This may react with copper or lead plumbing to form explosive metal azides. Upon disposal, flush with large amounts of water to prevent azide build up.
- Do not use the reagent after the expiration date printed on the kit.

H315 Causes skin irritation

H319 Causes serious eye irritation

P302+352: IF ON SKIN: Wash with soap and water

P305+P351+P338 IF IN EYES: Rinse continuously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing.

Wear suitable protective equipments

REAGENT PREPARATION

Reagent is supplied ready to use.

REAGENT STORAGE

Store reagent at room temperature. (2-25°C)

REAGENT DETERIORATION

The reagent should not be used if:

- The reagent is turbid.
- The reagent fails to meet stated parameters of performance

SPECIMEN COLLECTION AND HANDLING

Collection: No special preparation of the patient is necessary and sample preservatives are not required.

Sample Type: Serum, is the recommended specimen.

Sample Storage: Albumin in serum is stable for one month at 2-8 °C.

INTERFERENCES

Bilirubin: No interference up to 40 mg/dL.

Hemoglobin: No interference up to 400 mg/dL.

Lipemia: No interference in the presence of triglycerides up to 1500 mg/dL.

Ascorbic acid: No interference up to 25 mg/dL.

ADDITIONAL EQUIPMENT REQUIRED BUT NOT PROVIDED:

- A clinical chemistry analyzer capable maintaining constant temperature (37°C), and measuring absorbance at 630nm.
- Iron-free deionized water and related equipment, e.g.: pipettes
- Analyzer specific consumables, e.g.: sample and read cups
- Control and calibrator materials

PROCEDURE

Wavelength	: 630 nm
Working temperature:	: 37°C
Optical path	: 1 cm
Assay type	: Endpoint
Direction	: Increasing

	Blank	Standard	Sample
Reagent	1000 µL	1000 µL	1000 µL
Distilled water	10 µL	--	--
Standard	--	10 µL	--
Sample	--	--	10µL

Mix and then incubate for 90 second at 37°C. Measure the absorbance of sample and standard against the reagent blank.

Procedure Notes

- The reagent and sample volumes may be altered proportionally to accommodate different instruments requirements.
- The temperature of the reaction is not critical, however the temperature should be held constant.
- Unit Conversion : g/dL x 10 = g/L

CALCULATION

(Abs = Absorbance)

$\frac{\text{Abs. of Sample}}{\text{Abs. of standard}} \times \text{conc. of std. (g/dL)} = \text{ALBUMIN (g/dL)}$

Abs. of standard

LIMITATIONS

Samples with values exceeding 6.0 g/dL should be diluted 1:1 with saline and re-run. The final answer should be multiplied by two.

CALIBRATION

Use an aqueous albumin standard or an appropriate serum calibrator.

QUALITY CONTROL

The integrity of the reaction should be monitored by use of a two level control with known albumin values.

EXPECTED VALUES

3.5 – 5.0 g/dL

It is strongly recommended that each laboratory establish its own normal range.

PERFORMANCE

Sensitivity/ Limit of Detection(LOD):

The lower limit of detection is 0,025 g/dL.

Linearity:

When the recommendation is designed and studied, the test is linear up to 6 g/dL. Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

Precision :

	Within run (n=20)		Between run (n=20)	
Mean (g/dL)	3,16	5,26	3,29	5,36
SD	0,06	0,09	0,07	0,13
CV(%)	2,02	1,65	2,12	2,34

Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x).The results obtained using 40 samples were the following:

Correlation Coefficient (r) =0,9956

Regression $y=0,937x + 0,0659$

The results of the performance characteristics depend on the analyzer used.

REFERENCES

- Tietz, N., Textbook of Clinical Chemistry, Philadelphia, W.B. Saunders, 1986 pp. 701-704
- Webster D: 177. The Immediate Reaction between Bromocresol Green and Serum as a Measure of Albumin Content. Clin Chem 23:663
- Doumas BT, Warson WA, and Biggs NG: 1971. Albumin Standards and the measurement of Serum Albumin with Bromocresol Green Clin Chem Acta 31:87
- Young DS, Effects of Drugs on Clinical Laboratory Tests. Third Edition 1990;12-6

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BİLİMSEL TIBBİ ÜRÜNLER PAZ.SAN.VE TİC.LTD.ŞTİ