ALP-10100
ALP-10250
ALP-10500
ALP-10160A

2x40mL+2x10mL
5x40mL+1x50mL
5x80mL+1x100mL
4x32mL+4x8mL

ALP-10228M
ALP-10228P
ALP-10300M2
ALP-10200M3

4x45mL+2x24mL
4x45mL+2x24mL
6x40mL+3x20mL
4x40mL+2x20mL

ALP-10100A2

4x20mL+2x10mL

2-8 °C

KULLANIM AMACI

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BİLİMSEL KAYNAKÇA

- Wenger C. et al. Alkaline phosphatase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1094-1098.
- Rosalki S et al. Clin Chem 1993; 39/4: 648-652.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACCC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACCC Press, 1995.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACCC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACCC 1995.

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2-8 °C

INTENDED USE

For the in vitro quantitative determination of Alkaline Phosphatase in human serum.

METHODOLOGY

Alkaline phosphatase in serum is determined by measuring the rate of hydrolysis of various phosphate esters under specified conditions.

Principle:

Alkaline phosphatase (ALP) catalyses the hydrolysis of p-nitrophenyl phosphate at suitable pH, liberating p-nitrophenol and phosphate, according to the following reaction:

The rate of p-Nitrophenol formation, measured photometrically, is proportional to the catalytic concentration of alkaline phosphatase present in the sample.

REAGENT COMPOSITION

REAGENT 1

Diethanolamine 1.0 mM
Magnesium chloride 0.5 mM

REAGENT 2

p-Nitrophenylphosphate 10.0 mM
pH 9,7 ± 0.2

Precautions and Warnings:

- For in vitro diagnostic use only.
- DO NOT pipette by mouth. Avoid contact with skin and eyes. If spilt, thoroughly wash affected area with water.
- Reagent contains Sodium Azide as a preservative. In a dry state may react with copper or lead plumbing to form explosive metal azides. Upon disposal, flush with large amounts of water to prevent azide build up.
- Do not use the reagent after the expiration date printed on the kit.

REAGENT 1:

H315 Causes skin irritation

H319 Causes serious eye irritation

H373 May damage to organs through prolonged or repeated exposure.

P302+352: IF ON SKIN: Wash with soap

and water

P305+P351+P338 IF IN EYES: Rinse continuously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing

Wear suitable protective equipments

REAGENT PREPARATION

Reagents are ready to use, liquid form.

Working Reagent: The reagents can be combined to make a working solution by mixing 4 parts of Reagent 1 with 1 part of Reagent 2 (e.g., 20mL Rgt 1 to 5 mL Rgt 2).

STABILITY AND STORAGE

- Store the reagents at 2-8°C (refrigerated).
- The reagents are stable until the expiration date when stored at 2-8°C.
- Working reagent is stable for 4 weeks when stored tightly closed at (2-8°C), or 10 days at room temperature protected from light and contamination during their use.
- Do not freeze the reagents.
- Reagent 2 should be protected from light.

REAGENT DETERIORATION

Do not use the reagent if:

- The reagent is turbid and presence particles.
- The reagent has an optical density greater than 1.3 at 405 nm.

SPECIMEN COLLECTION AND STORAGE

- Use non-hemolyzed serum or heparin plasma
- Serum samples should be stored 3 days at 2-8 °C.

INTERFERENCES

Fluoride, oxalate, citrate and EDTA inhibit alkaline phosphate activity and should therefore not be used as anticoagulants.

Bilirubin: No interference up to 40 mg/dL.

Hemoglobin: No interference up to 100 mg/dL.

Lipemia: No interference in the presence of triglycerides up to 1500 mg/dL.

Ascorbic acid: No interference up to 25 mg/dL.

ADDITIONAL EQUIPMENT REQUIRED BUT NOT PROVIDED

- A clinical chemistry analyzer capable maintaining constant temperature (37°C), and measuring absorbance at 405nm.
- Deionized water and related equipment, e.g.: pipettes

3. Analyzer specific consumables, e.g.: read and sample cups

4. Control materials.

PROCEDURE

Wavelength : 405 nm
Working temperature : 37°C
Optical path : 1cm
Assay type : Kinetic
Direction : Increasing

1. Bring reagent to room temperature (15-30°C).

2. Set the photometer to 0 (zero) absorbance with distilled water.

	Sample
Reagent 1	800 µL
Reagent 2	200 µL
Sample	20 µL

3. After 2-3 minutes, read and record absorbance. Return tube to 37°C.

4. Repeat readings every minute for the next three minutes.

5. Calculate the average absorbance difference per minute. ((Abs/Min.)

CALCULATIONS

One Unit (U/L) is defined as the amount of enzyme that catalyzes the transformation of one micromole of substrate per minute under the specified conditions.

ALP U/L = Δ Abs./min. x 3300

Temperature conversion factors

To correct results to other temperatures multiply by:

Assay temperatures	Conversion factors to		
	25°C	30°C	37°C
25°C	1.00	1.22	2.64
30°C	0.82	1.00	1.33
37°C	0.61	0.75	1.00

QUALITY CONTROL:

The integrity of the reagent should be monitored by use of a two level control with known Alkaline Phosphates values.

EXPECTED VALUES

	25°C	30°C	37°C
Children (1-14 years)	< 400 U/L	< 480 U/L	< 645 U/L
Adults	60-170 U/L	73-207 U/L	98-279 U/L

It is strongly suggested that each laboratory establish its own normal range.

PERFORMANCE

Sensitivity/ Limit of Detection(LOD):

The lower limit of detection is 2,63 U/L

Linearity:

When the recommendation is designed and studied, the test is linear up to 1000 U/L. Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

Precision :

	Within run (n=20)		Between run (n=20)	
Mean (U/L)	94,4	227,9	94,1	226,9
SD	1,13	2,42	0,95	1,30
CV(%)	1,20	1,06	1,01	0,57

Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x). The results obtained using 40 samples were the following:

Correlation Coefficient (r) =0,9985

Regression y =0,9026 x +16,473

The results of the performance characteristics depend on the analyzer used.

REFERENCES

- Wenger C. et al. Alkaline phosphatase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1094-1098.
- Rosalki S et al. Clin Chem 1993; 39/4: 648-652.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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BİLİMSEL TIBBİ ÜRÜNLER PAZ.SAN.VE TİC.LTD.ŞTİ