REF	CONT	REF	CONT	REF	CONT
AMY-10050	1x40mL+1x10mL	AMY-10160A	4x32mL+4x8 mL	AMY-10200M3	4x40mL+2x20mL
AMY-10100	2x40mL+2x10mL	AMY-10152M	4x30mL+2x16mL	AMY-10040A2	4x8mL+1x8 mL
AMY-10250	5x40mL+1x50 mL	AMY-10152P	4x30mL+2x16mL		
AMY-10500	5x80mL+1x100mL	AMY-10100M2	4x20mL+2x10mL		

: 2-8 °C

**KULLANIM AMACI**

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BİLİMSEL KAYNAKÇA

1. Tietz, N. W., et al. Abs. of Proc. Of Int'l Seminar and workshop on Enzymology, Chicago, IL (May 972).
2. Morishita, et al. Clinical Chemistry 46 : 7, 928 - 933.
3. Young, DS., Effects of Drugs on Clinical Laboratory Tests, fifth edition 2000, AACCC Press, Washington, D.C.
4. Tietz , N.W. clinical Guide to laboratory tests., Philadelphia , W.B.Saunders company ,p.54 (1983)
5. Tietz, N. W. Textbook of clinical chemistry, Philadelphia, W.B. Saunders Com-pany, pp. 725 - 734 (1986).
6. Tietz , N. W., Fundamentals of Clinical Chemistry, Philadelphia, W.B. Saunders Company, p. 54 (1983).IFCC methods for the measurement of catalytic concentrations of enzymes, J. Clin. Chem. Acta. (1999) 281: 5.
7. Clinical Chemistry, Principles, Procedures, Correlations, Michael L.
8. Bishop et. al., 5th Edition.
9. Data on file: Coral Clinical Systems

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

INTENDED USE

For the quantitative kinetic determination of amylase activity in human serum using a manual or automated procedure.

Principle:

α -Amylase hydrolyzes the substrate (Gal-G2- α -CNP) 2-chloro-4-nitrophenyl-L-O- α -D-galactopyranosyl maltoside to release directly 2-chloro -p-nitrophenol (CNP), which is measured colorimetrically at 405 nm. The rate of increase in absorbance is proportional to α -amylase activity in the sample.

REAGENT COMPOSITION

Reagent R1

MES Buffer	50	mmol/L
Calcium chloride	6	mmol/L
Sodium chloride	300	mmol/L
KSCN	140	mmol/L
preservative and stabilizer		
pH = 6.0		

Reagent R2

MES Buffer	50	mmol/L
Gal-G2- α -CNP	4.0	mmol/L
pH = 6.0		
Preservatives		

PREPARATION OF WORKING REAGENT

R1: Liquid, R2: Liquid

Prepare working reagent by mixing 4 parts reagent 1 with 1 part reagent 2 (e.g. 800 μ L R1 with 200 μ L R2 reagent)

REAGENT STORAGE

1. Store R1 and R2 reagents at 2-8°C.
2. Store working reagent at 2-8°C.
3. Working reagent is stable 2 weeks at 2-8°C.
4. Do not freeze the reagents

REAGENT DETERIORATION

Do not use if:

1. The absorbance of the working reagent is greater than 0.600 when measured at 405 nm against water in a cuvette with a 1 cm path length.
2. The reagent fails to meet stated parameters of performance.
3. The reagent is turbid or displays other evidence of bacterial contamination.

PRECAUTIONS

1. This reagent kit is intended for *in vitro* diagnostic use only.
2. This reagent contains potassium thiocyanate. Avoid ingestion.
3. Reagents are acidic solutions, flush with water when contact occurs.
4. This reagent contains sodium azide (0.01%) as preservative. Do not ingest. May react with lead and copper plumbing to form highly explosive metal azides. Upon disposal, flush with a large volume of water to prevent azide build up.

Reagent 1 and Reagent 2 contains Sodium hydroxide .

R36/38 : Irritation to eyes and skin.

In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective cloth gloves and eye/face protection.

SPECIMEN COLLECTION AND HANDLING

1. Unhemolyzed serum is the specimen of choice.
2. Anticoagulants, such as Citrate and EDTA, bind calcium that is needed for amylase activity. Plasma with these anticoagulants should not be used.
3. Amylase in serum is reported stable for one week at room temperature (18-25 °C) and for one month when stored refrigerated at 2-8 °C.

INTERFERENCES

A number of drugs and substances affect the determination of amylase.

Bilirubin: No interference up to 40 mg/dL.

Hemoglobin: No interference up to 500 mg/dL.

Lipemia: No interference in the presence of triglycerides up to 1500 mg/dL.

Ascorbic acid: No interference up to 25 mg/dL.

MATERIALS REQUIRED BUT NOT PROVIDED

1. Accurate pipetting devices
2. Test tubes/rack/Timer
3. Heating block or water bath (37°C).
4. Spectrophotometer capable to read at 405nm. The cuvette should be temperature controlled to maintain temperature (37°C) during the assay.

PROCEDURE

Wavelength	: 405 nm
Working temperature	: 37°C
Optical path	: 1cm
Assay type	: Rate/kinetic
Direction	: Increasing
Bring reagent to room temperature (15-30°C).	

1. Set the photometer to 0 (zero) absorbance with distilled water.

Reagent 1	Sample
800 μ L	800 μ L
200 μ L	200 μ L
25 μ L	25 μ L

1. After one minute, read and record absorbance. Return tube to 37°C.
2. Repeat readings every minute for the next two minutes
3. Calculate the average absorbance difference per minute. (Δ Abs/Min.)
4. The Δ Abs/Min multiplied by the factor 3178 (See Calculation) will yield results in IU/L.

CALCULATION

Δ Abs./min x 3178 = U/L α -amylase

CALIBRATION

For the calibration of automated photometric systems the bt products Chem-Calibrator is recommended.

QUALITY CONTROL

Internal quality control we recommend bt products Chem-Control Normal and Chem-Control Abnormal controls to be assayed.

EXPECTED VALUES

Serum: 25-125 U/L

Since the expected values are affected by age, sex, diet and geographical location, each laboratory is strongly urged to establish its own reference range for this procedure.

PERFORMANCE

Sensitivity/ Limit of Detection(LOD):

The lower limit of detection is 0,596 U/L

Linearity:

When the recommendation is designed and studied, the test is linear up to 1800 U/L. Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

Precision :

	Within run (n=20)		Between run (n=20)	
Mean (U/L)	80,1	178,0	79,0	183,3
SD	0,72	1,80	0,87	2,70
CV(%)	0,90	1,01	1,10	1,47

Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x). The results obtained using 40 samples were the following:

Correlation Coefficient (r) =0,9998

Regression $y=1,001x + 5,5164$

The results of the performance characteristics depend on the analyzer used

REFERENCES

1. Tietz, N. W., et. al. Abs. of Proc. Of Int'l Seminar and workshop on Enzymology, Chicago, IL (May 972).
2. Morishita, et al. Clinical Chemistry 46 : 7, 928 - 933.
3. Young, DS., Effects of Drugs on Clinical Laboratory Tests, fifth edition 2000, AACCPress, Washington, D.C.
4. Tietz, N.W. clinical Guide to laboratory tests.,
5. Philadelphia, W.B.Saunders company, p.54 (1983)
6. Tietz, N. W. Textbook of clinical chemistry, Philadelphia, W.B. Saunders Com-pany, pp. 725 - 734 (1986).
7. Tietz, N. W., Fundamentals of Clinical Chemistry, Philadelphia, W.B. Saunders Company, p. 54 (1983).IFCC methods for the measurement of catalytic concentrations of enzymes, J. Clin. Chem. Acta. (1999) 281: 5.
8. Clinical Chemistry, Principles, Procedures, Correlations, Michael L.
9. Bishop et. al., 5th Edition.
10. Data on file: Coral Clinical Systems

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BİLİMSEL TIBBİ ÜRÜNLER PAZ.SAN.VE TİC.LTD.ŞT