REF	CONT	REF	CONT
ASO-30050	1x40mL+1x10mL	ASO-30150M	4x30mL+2x15mL
ASO-30150	3x40mL+1x30mL	ASO-30150P	4x30mL+2x15mL
SO-30120A	4 x24mL+4x6mL	ASO-30100M2	2x40mL+2x10mL
ASO-30050A2	4x10mL+1x10mL	ASO-30200M3	4x40mL+2x20mL

: 2-8 °C

**KULLANIM AMACI**

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
**INTENDED USE**

Quantitative determination of anti-streptolysin O (ASO).

CLINICAL SIGNIFICANCE

SLO is a toxic immunogenic exoenzyme produced by β -hemolytic Streptococci of groups A, C and G. Measuring the ASO antibodies are useful for the diagnostic of rheumatoid fever, acute glomerulonephritis and streptococcal infections. Rheumatic fever is an inflammatory disease affecting connective tissue from several parts of human body as skin, heart, joints etc... and acute glomerulonephritis is a renal infection that affects mainly to renal glomerulus.

Principle:

The ASO-Turbilatex is a quantitative turbidimetric test for the measurement of ASO in human serum or plasma.

Latex particles coated with streptolysin O (SLO) are agglutinated when mixed with samples containing ASO. The agglutination causes an absorbance change, dependent upon the ASO contents of the patient sample that can be quantified by comparison from a calibrator of known ASO concentration.

REAGENTS

Diluent (R1)	Tris buffer 20 mmol/L, pH 8.2. Sodium azide 0.90 g/L. Merthiolate 0.05 g/L
Latex (R2)	Latex particles coated with streptolysin O, pH 10.0. Sodium azide 0.90 g/L, Merthiolate 0.05 g/L
Optional	ASO Calibrator. ASO/CRP/RF Control Level 1 ASO/CRP/RF Control Level 2

PRECAUTIONS

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

bt products recommends using ASOCAL-30 calibrator .

Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

REAGENT PREPARATION

Reagents are supplied in a two vial, ready to use, liquid form.

Swirl the latex vial gently before use.

For some analyzers, the reagents can be combined to make a working solution by mixing 4 part of Diluent (R1) with 1 part of Latex (R2) (e.g., 4mL R1 + 1mL R2).

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations prevented during their use. Do not use reagents over the expiration date. Do not freeze

Working reagent: Stable for 30 days at 2-8°C.

Reagent deterioration:

Presence of particles and turbidity.

Frozen Latex or Diluent could change the functionality of the test.

INTERFERENCES

Bilirubin (20 mg/dL), hemoglobin (10 g/L), lipemia (10 g/L) and rheumatoid factors (600 IU/mL), do not interfere. Other substances may interfere.

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer thermostatable at 37°C with a 540(530-550)nm filter.

SAMPLES

Fresh serum. Stable 7 days at 2-8°C or 3 months at (-20°C).

Samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolized or lipemic samples.

PROCEDURE

Wavelength	: 540(530-550)nm
Working temperature	: 37°C
Optical path	: 1 cm
Assay type	: Turbidimetry
Direction	: Increasing

1. Bring the reagents and the photometer (cuvette holder) to 37°C .
2. Adjust the instrument to zero with distilled water.
3. Pipette into a test tube

Sample or calibrator tube	
R1 (DILUENT)	800 μ L
R2 (LATEX)	200 μ L
Sample or calibrator	10 μ L

4. Mix well and immediately read the absorbance (Abs.1). Read the absorbance (Abs.2) exactly 2 minutes after the sample or calibrator addition.

CALCULATIONS

$(\text{Abs.2}-\text{Abs.1}) \text{ sample} \times \text{Calibrator concentration} = \text{ASO IU/mL}$
 $(\text{Abs.2}-\text{Abs.1}) \text{ calibrator}$

QUALITY CONTROL

ASO/CRP/RF Control Level 1 or ASO/CRP/RF Control Level 2 Control sera are recommended to monitor the performance of manual and automated assay procedures. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES

Adults : < 200 IU/mL

Children < 5 years : < 100 IU/mL

Each laboratory should establish its own reference range.

PERFORMANCE

1. **Linearity limit:** Up to 800 IU/mL, under the described assay conditions. Samples with higher concentrations, should be diluted 1/3 in NaCl 0.9 g/dL and retested again. The linearity limit depends on the sample-reagent ratio, as well the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
2. **Detection limit:** Values less than 20 IU/mL give non-reproducible results.
3. **Prozone effect:** No prozone effect was detected up to 3000 IU/mL.
4. **Sensitivity:** Δ 0.73 mA. IU/mL.
5. **Precision:**

EP5(20 day)	% CV
+/- 100 IU/mL	
+/- 200 IU/mL	
+/- 400 IU/mL	
Total	6.4 % 5.7 % 5.1 %
Within Run	2.4 % 1.7 % 1.4 %
Between	3.6 % 4.2 % 4.9 %
Run	
Between	4.7 % 3.5 % 0.7 %
Day	
6. **Accuracy:** Results obtained using this reagent (y) were compared to those obtained using a commercial reagent (x) with similar characteristics. 60 samples ranging from 20 to 800 IU/mL of ASO were assayed. The correlation coefficient (r) was 0.99 and the regression equation $y = 0.915x - 4.844$.

The results of the performance characteristics depend on the analyzer used.

NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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Bilimsel Tıbbi Ürünler Paz.San.Tic.Ltd.Şti.