REF	CONT	REF	CONT	REF	CONT
BIT-50100	2x40 mL+ 2x10 mL	BIT-50228M	4x45 mL +2x24 mL		
BIT-50250	5x40 mL+1x50 mL	BIT-50228P	4x45 mL+2x24 mL		
BIT-50500	5x80 mL+1x100 mL	BIT-50200M2	4x40 mL+2x20 mL		
BIT-50160A	4x32 mL +4x8 mL	BIT-50240M3	4x48 mL+2x24 mL		

⌘ : 2-8 °C

KULLANIM AMACI

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⌘ : 2-8 °C


INTENDED USE

For the quantitative *in vitro* determination of total bilirubin in serum.

PRINCIPLE OF THE METHOD

When a sample is mixed with the reagent containing the detergent and the vanadate, at around pH 3.0, total bilirubin in the sample is oxidized to biliverdin. This causes the absorbance of yellow, specific to bilirubin, to decrease. Therefore, the total bilirubin concentration in the sample can be obtained by measuring the absorbances before and after vanadate oxidation.

REAGENT COMPOSITION

Active Ingredients

REAGENT 1 (Buffer Solution)

Citrate Buffer pH 2.9 0.1 mol/L
Detergent

REAGENT 2 (Vanadate Solution)

Phosphate Buffer pH 7,0 10 mmol/L
Sodium Metavanadate 4 mmol/L

Precautions:

1. This reagent is for *in vitro* diagnostic use only.
2. Not to be used internally in humans or animals.
3. Do not mix the reagents from one test unit with those of another test unit which has different lot number.
4. Do not use reagents past the expiration date stated on each reagent container label.
5. Do not use the reagents for any other purpose than described herein.
6. After opening the reagent, it is not recommended to store it for a long period of time. When the opened reagent is stored, cap the bottle and keep it at the specified temperature.
7. When discarding the reagents, dispose of them according to local or national regulations.

REAGENT 1:

H319 Causes serious eye irritation

H351 Suspected to causing cancer

P308+P313: If exposed or concerned, get medical advice/attention

P305+P351+P338 IF IN EYES :R Inse continuously with water for several minutes .Remove contact lenses if present and easy to do. Continue rinsing

Wear suitable protective cloth

Use appropriate container to avoid environmental contamination

REAGENT PREPARATION

Reagent is supplied ready to use.

REAGENT STORAGE

Packaged reagents are stable until expiration date when stored at 2-8°C.

REAGENT DETERIORATION

Do not use the reagent if:

1. The presence of precipitates in the reagents, turbidity.
2. If the values of control sera are outside the manufacturer's acceptable range.

SPECIMEN COLLECTION AND STORAGE

Fresh, unhemolyzed serum is recommended. When stored the serum must be frozen (-20°C) under conditions with no light exposure since serum bilirubin degrades to biliverdin by light.

INTERFERENCES

Lipemia: No interference in the presence of triglycerides up to 1500 mg/dL.

Ascorbic acid: No interference up to 15 mg/dL.

ADDITIONAL EQUIPMENT REQUIRED BUT NOT PROVIDED

1. A clinical chemistry analyzer capable maintaining constant temperature (37°C), and measuring absorbance at (450nm - 546 nm).
2. Deionized water and related equipment, e.g.: pipettes
3. Analyzer specific consumables, e.g.: sample cups
4. Control and calibrator materials.

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Rev.Date / No: 03.10.2022/6

	Blank	Standard	Sample
Reagent 1	280 µL	280 µL	280 µL
Sample	--	-	10 µL
Standard	--	10 µL	--
Distilled Water	10 µL	--	--

Mix and then incubate for 5 minutes at 37°C. Read Abs.1.Then add Reagent 2 to all tubes.

Reagent 2	70 µL	70 µL	70 µL
-----------	-------	-------	-------

Mix and incubate 5 minutes at 37°C. Read Abs.2

Procedure Notes:

1. The reagent and sample volumes may be altered proportionally to accommodate various instrument requirements.
2. Bilirubin is extremely light sensitive. All samples should be stored protected from light sources.

CALCULATIONS:

(Abs.= Absorbance)

$$\frac{\text{Sample (Abs.2 - Abs.1)}}{\text{Standard (Abs.2 - Abs.1)}} \times \text{Conc. of Standard} = \text{Total Bilirubin (mg/dL)}$$

LIMITATIONS:

Samples with values exceeding 40 mg/dL should be diluted 1:1 with saline and re-run. The final answer should be multiplied by two.

CALIBRATION

Use a **bt products** Total Bilirubin standard or an appropriate serum calibrator.

QUALITY CONTROL

The integrity of the reaction should be monitored by use of a two level control with known Total Bilirubin values.

EXPECTED VALUES

Adults and infants over 1 month old : 0.2 – 1.0 mg/dL.

It is recommended each laboratory verify this range or derives a reference interval for the population that it serves.

PERFORMANCE

Sensitivity/ Limit of Detection(LOD):

The lower limit of detection is 0,066 mg/dL.

Linearity:

When the recommendation is designed and studied, the test is linear up to 30 mg/dL. Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

Precision :

	Within run (n=20)	Between run (n=20)		
Mean (U/L)	1,00	4,54	1,01	4,71
SD	0,04	0,04	0,03	0,03
CV(%)	4,35	0,82	2,49	0,68

Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x).The results obtained using 40 samples were the following:

Correlation Coefficient (r) =0,9967

Regression $y=0,9362x - 0,0542$

The results of the performance characteristics depend on the analyzer used.

REFERENCES

1. Malloy H.T.Evelyn K.L.The determination of bilirubin with the photoelectric colorimetry. J.Jbiol.Chem.199:481-490,(1937).
2. Jendrassik L.Cleghom R.A.Photometrische bilirubinbestimmung.Biochem.Z.289:1-14,(1937)
3. Michaelson M.Bilirubin determination in serum and urine.Scand J.Clin.Lab.Invest.12(Suppl 56):1-80,(1937).
4. Murao S.Tanaka N.A new enzyme bilirubin oxidase,produced by Myrothecium varrucaria MT-1.Agric.Biol.Chem.45:2383-2384,(1981)
5. Tokuda K.and tanimoto K.New method of measuring serum bilirubin using vanadic acid.Jpn.J.Clin.Chem.22(2),116-122,(1993)
6. Carl A.Burtis,Edward R.Ashwood,Tietz text book of clinical chemistry,second edition,P1468(SAUNDERS)