REF	CONT	REF	CONT	REF	CONT
CAL-10100	2 x 50 mL	CAL-10180M	4 x 45 mL	CAL-10100A2	5 x 20 mL
CAL-10300	6 x 50 mL	CAL-10180P	4 x 45 mL		
CAL-10600	6 x 100 mL	CAL-10240M2	6 x 40 mL		
CAL-10080A	4 x 20 mL	CAL-10240M3	6 x 40 mL		

: 15-25 °C

## KULLANIM AMACI

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

### INTENDED USE

For the *in vitro* quantitative determination of Calcium in serum.

### METHODOLOGY

This method uses Arsenazo III and has been modified to provide a highly sensitive and stable reagent system

### Principle:

Calcium reacts with Arsenazo III in a slightly alkaline medium to form a purple-colored complex which absorbs at 650 nm. The intensity of the color is proportional to the calcium concentration.

### REAGENT COMPOSITION

Active Ingredients	Concentration
Arsenazo III	0.2 mM
Imidazol Buffer	100 mM
Surfactant	
pH 8.5	

### Precautions

1. Reagent is for "in vitro" diagnostic use only.
2. Reagent may be irritating to the skin. Avoid contact. Flush with water if contact occurs.
3. Reagent contains Sodium Azide as a preservative.  
*H360d May damage fertility or the unborn child*  
*P308+313 If exposed or concerned: get medical advice / attention*

### REAGENT PREPARATION

Reagent is supplied ready to use.

### STABILITY AND STORAGE

When stored at 15-25°C, the reagent is stable until the expiration date stated on the label.

### REAGENT DETERIORATION

The reagent should not be used if:

1. The reagent is turbid.
2. The reagent fails to meet stated parameters of performance.

### SPECIMEN COLLECTION AND STORAGE

1. Fresh, unhemolyzed serum is the preferred specimen.
2. Anticoagulants other than heparin should not be used.
3. Remove serum from clot as soon as possible since red cells can absorb calcium.
4. Older serum specimens containing visible precipitate should not be used.
5. Serum calcium is stable for one week at refrigerated (2-8°C), and up to five months frozen (-15 to -25 °C) when protected from evaporation. Specimens should not be thawed and then refrozen.

### INTERFERENCES

Bilirubin: No interference up to 40 mg/dL.

Hemoglobin: No interference up to 500 mg/dL.

Lipemia: No interference in the presence of triglycerides up to 1500 mg/dL.

Ascorbic acid: No interference up to 25 mg/dL.

### ADDITIONAL EQUIPMENT REQUIRED BUT NOT PROVIDED

1. A clinical chemistry analyzer capable maintaining constant temperature (37°C), and measuring absorbance at 650 nm (600 – 660).
2. Deionized water and related equipment, e.g.: pipettes
3. Analyzer specific consumables, e.g.: sample and read cups
4. Control and calibrator materials.

### PROCEDURE

Wavelength:	:650 nm
Working temperature:	:37°C
Assay type:	:Endpoint
Direction:	:Increasing

	Blank	Standard	Sample
Reagent	1000 µL	1000 µL	1000 µL
Distilled water	10 µL	--	--
Standard	--	10 µL	--
Sample	--	--	10µL

Mix and then incubate for min 60 second at 37°C. Measure the absorbance of sample and standart against the reagent blank.

### Procedure Notes

1. Final color is stable for 60 minutes.
2. Samples with calcium above 15 mg/dL should be diluted 1:1 with saline, re-assayed, and the result multiplied by two.
3. Severely lipemic samples may require a serum blank.
4. Contamination of glassware with calcium will adversely affect test results. Acid-washed glass or plastic test tubes are recommended.
5. Reagent and sample volumes may need to be adjusted according to the requirements of the specific spectrophotometer being used.

### CALCULATION

(Abs = Absorbance)

$\text{Abs. Sample} \times \text{Conc. of Standard} = \text{Calcium (mg/dL)}$

Abs. Standard

### CALIBRATION

Use an aqueous Calcium standard or an appropriate serum calibrator.

### QUALITY CONTROL

The integrity of the reagent should be monitored by use of a two level control with known Calcium values.

### EXPECTED VALUES

Adults: 8.5 - 10.4 mg/dL

Newborns: 7.8 - 11.2 mg/dL

It is strongly recommended that each laboratory establish its own normal range.

### PERFORMANCE

#### Sensitivity/ Limit of Detection(LOD):

The lower limit of detection is 0,99 mg/dL.

#### Linearity:

When the recommendation is designed and studied, the test is linear up to 17 mg/dL. Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

#### Precision :

	Within run (n=20)		Between run (n=20)	
Mean (mg/dL)	9,4	14,08	9,19	13,29
SD	0,15	0,21	0,20	0,15
CV(%)	1,55	1,51	2,21	1,10

#### Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x).The results obtained using 40 samples were the following:

Correlation Coefficient (r) =0,9968

Regression  $y=1,0314x - 0,0291$

The results of the performance characteristics depend on the analyzer used.

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Bilimsel Tıbbi Ürünler Paz.San.Tic.Ltd.Şti.