REF	CONT	REF	CONT
CHE-10100	2x42 mL+2x8,5 mL	CHE-10050M2	2x21 mL+2x4,5 mL
CHE-10050A	2x21 mL+2x4,5 mL	CHE-10100M3	4x21 mL+2x8,5 mL
CHE-10075M	3x21mL+2x6,5 mL	CHE-10060A2	4x12,5 mL+2x5 mL
CHE-10075P	3x21mL+2x6,5 mL		

: 2-8 °C

**KULLANIM AMACI**

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
**INTENDED USE**



For the in vitro quantitative determination of Cholinesterase in serum or plasma.

**METHODOLOGY**

Kinetic photometric test, optimized method according to the recommendation of the German Society of Clinical Chemistry (DGKC).

**Principle:**

Cholinesterase hydrolyses butyrylthiocholine under release of butyric acid and thiocholine. Thiocholine reduces yellow potassium hexacyanoferrate (III) to colorless potassium hexacyanoferrate (II). The decrease of absorbance is measured at 405 nm.

**REAGENT COMPOSITION**

<b>Reagent 1</b>	:Pyrophosphate buffer pH 7,6	75 mmol/L
	Potasyum hexacyanoferrate (III)	2 mmol/L
<b>Reagent 2</b>	:Butyrylthiocholine	>12 mmol/L

**REAGENT PREPARATION**

Ready to use

**REAGENT STORAGE & STABILITY**

The reagents are stable up to the end of the indicated month of expiry, if stored at 2 – 8 °C and contamination is avoided. Do not freeze the reagents and protect from light.

**WARNING AND PRECAUTIONS**

1. In very rare cases, samples of patients with gammopathy might give falsified results.
2. Please refer to the safety data sheets and take the necessary precautions for the use of laboratory reagents. For diagnostic purposes, the results should always be assessed with the patient's medical history, clinical examinations and other findings.
3. For professional use only.

**SPECIMEN COLLECTION AND STORAGE**

1. Serum, plasma with heparine or EDTA.
2. Stability 1 weeks at 2 – 8 °C , 1 year at – 20 °C

**INTERFERENCES**

Bilirubin: No interference up to 40 mg/dL.  
Hemoglobin: No interference up to 500 mg/dL.  
Lipemia: No interference in the presence of triglycerides up to 1500 mg/dL.  
Ascorbic acid: No interference up to 25 mg/dL.

**ADDITIONAL EQUIPMENT REQUIRED BUT NOT PROVIDED**

1. A clinical chemistry analyzer capable maintaining constant temperature 37 °C, and measuring absorbance at 405 nm.
2. Deionized water and related equipment, e.g.: pipettes
3. Analyzer specific consumables, e.g.: sample cups
4. Control and calibrator materials.
5. NaCl ( 9 g/L) solution

**PROCEDURE**

Wavelength	: 405 nm
Temperature	: 37 °C
Optical path	: 1cm
Assay type	: Against reagent blank

Sample/calibrator	Reagent blank	Sample/calibrator
Distilled water	--	4 µL
Reagent 1	4 µL	--
Mix and incubate 5 min.at 37°C.Read absorbance then add.	200 µL	200 µL
Reagent 2	40 µL	40 µL
Mix and incubate for 90 sec. at 37°C.The absorbance within 100seconds was continuously measured to calculate the change rate of absorbance (ΔA/min).		

**CALCULATIONS**

$$\frac{\Delta A / \text{min.sample}}{\Delta A / \text{min calibrator}} \times \text{con.of calibrator(U/L)} = \text{CHE (U/L)}$$

**CONVERSION FACTOR**

$$\text{CHE (KU/L)} \times 16.67 = \text{CHE (µkat/L)}$$

**QUALITY CONTROL**

The integrity of the reaction should be monitored by use of a two level control with known Cholinesterase values.

**EXPECTED VALUES**

Men : 5,1 – 11,7 KU/L (5.100-11.700 U/L)  
Women : 4,0 – 12,6 KU/L (4.000- 12.600 U/L)

It is highly recommended that each laboratory establish its own expected range.

**PERFORMANCE****Sensitivity/ Limit of Detection(LOD):**

The lower limit of detection is 0,068 U/L.

**Linearity:**

When the recommendation is designed and studied, the test is linear up to 20.000 U/L . Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

**Precision :**

	Within run (n=20)	Between run (n=20)		
Mean (U/L)	4787,7	5853,2	4787,1	5855,1
SD	8,50	7,05	9,26	8,49
CV(%)	0,18	0,12	0,19	0,14

**Accuracy (Method Comparison) :**

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x).The results obtained using 40 samples were the following:  
Correlation Coefficient (r) =0,9996  
Regression y=0,9992x +41,381  
The results of the performance characteristics depend on the analyzer used.

**REFERENCES**

- 1.Knedel, M., and R. Bottger. Klin. Wschr. 1967; 45: 325.
- 2.Han Zhijun, Huang Zhifeng, etc. Automatic analysis method for commonly used items in clinical chemistry. Liaoning: Liaoning Science and Technology Press, 2005. 1066-1103.

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BİLİMSEL TIBBİ ÜRÜNLER PAZ.SAN.VE TİC.LTD.ŞTİ