
KULLANIM AMACI

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
: 2-8 °C

**INTENDED USE**

For the quantitative determination of Chloride in human serum and plasma.

METHODOLOGY

The quantitative displacement of thiocyanate by chloride from mercuric thiocyanate and subsequent formation of a red ferric thiocyanate complex is measured colorimetrically:

The intensity of the color formed is proportional to the chloride ion concentration in the sample.

REAGENT COMPOSITION**Reagent**

Mercuric thiocyanate	4 mM
Ferric nitrate	40 mM
Mercuric nitrate	2 mM
Nitric acid	45 mM

Standard

SCHL-103	1x2 mL
----------	--------

REAGENT PREPARATION

All the reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

Precautions:*H314 Causes severe skin burns and serious eye damage.**P303+P361+P353 IF ON SKIN: Remove/take off immediately all contaminated clothing. Rinse skin with water/shower**P304+P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.**P305+P351+P338 IF IN EYES: Rinse continuously with water for several minutes**.Remove contact lenses if present and easy to do. Continue rinsing.**This material and its container must be disposed of as hazardous.**Avoid release to the environment.***REAGENT DETERIORATION**Presence of particles and turbidity. Blank absorbance (A) at 480 nm $\geq 0,15$ **SPECIMEN COLLECTION AND STORAGE**

- Sample: Serum, plasma, CSF or sweat. Sample without hemolysis and separated from cells as quickly as possible. Anticoagulants such as oxalate or EDTA are not acceptable they will interfere with results.
- Urine: Collect 24 hour urine specimen in chloride free containers. Dilute a sample 1/2 in distilled water. Mix. Multiply results by 2 (dilution factor).
- Stability of the sample : Ion chloride is stable 1 week at room temperature (15-25°C), in refrigerator (2-8°C), or frozen (-20°C) temperatures

INTERFERENCES**Bilirubin:** No interference up to 40 mg/dL.**Hemoglobin:** No interference up to 300 mg/dL.**Lipemia:** No interference in the presence of triglycerides up to 1500 mg/dL.**Ascorbic acid:** No interference up to 25 mg/dL.**Materials Required but not Provided**

- Spectrophotometer or colorimeter measuring at 480(440-500)nm.
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

PROCEDURE

Wavelength	: 480(440-500) nm
Assay Type	: Endpoint
Reaction Direction	: Increasing
Temperature	: 37°C or 15-25°C
Optical path	: 1 cm

- Bring to room temperature (15 -25 °C)
- Set the photometer to 0 (zero) with distilled water.

	BLANK	STANDARD	SAMPLE
Reagent	1000 µL	1000 µL	1000 µL
Standard	--	10 µL	--
Sample	--	--	10 µL

3. Mix and incubate for 5 min at 37°C/15-25°C

4. Read the absorbance (A) of the samples and standard, against the blank. The color is stable for at least 30 minutes.

CALCULATIONS
$$\frac{(\Delta A) \text{ Sample}}{(\Delta A) \text{ Standard}} \times \text{STD. conc.} = \text{mmol/L CHLORIDE in the sample}$$
Urine 24 h:
$$\frac{(\Delta A) \text{ Sample}}{(\Delta A) \text{ Standard}} \times \text{STD. conc.} \times \text{vol. (dL) urine / 24h} = \text{mmol/24h}$$

CHLORIDE in the Urine

Conversion factor: mmol/L = mEq/L

QUALITY CONTROL

The validity of the reaction should be monitored using control sera with known normal and abnormal values. If control values are found outside the defined range, check the instrument, reagent and calibration for problems.

Each laboratory should establish its own Quality control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES

Serum or plasma : 95 – 115 mmol/L

CSF : 95 – 110 mmol/L

Urine : 110 – 250 mmol/24h

Sweat : < 60 mmol/L

It is recommended that each laboratory establish its own reference range.

PERFORMANCE**Sensitivity/ Limit of Detection (LOD):**

The lower limit of detection is 1,70 mmol/L.

Linearity:

When the recommendation is designed and studied, the test is linear up to 160 mmol/L. Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

Precision :

	Within run (n=20)	Between run (n=20)		
Mean (mmol/L)	85,3	94,0	83,2	93,1
SD	1,16	1,19	1,18	1,04
CV(%)	1,36	1,27	1,41	1,12

Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x). The results obtained using 40 samples were the following:

Correlation Coefficient (r) = 0,9998

Regression $y = 1,001x + 5,5164$

The results of the performance characteristics depend on the analyzer used.

NOTES

- SCHL-103 CHLORIDE STD. : Proceed carefully with this product because due its nature it can get contaminated easily.
- It is recommended to use disposable material. If glassware is used the material should be scrupulously cleaned with H₂SO₄ - K₂Cr₂O₇ Solution and then thoroughly rinsed it with distilled water.
- Most of the detergents and water softening products used in the laboratories contains chelating agents. A defective rinsing will invalidate the procedure.
- Avoid the contact with metal materials.
- Use clean disposable pipette tips for its dispensation.

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www.btproducts.com.tr - www.bilimselitip.com

Rev. Date / No: 30.08.2023 / 6

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BİLİMSEL TIBBİ ÜRÜNLER PAZ.SAN.VE TİC.LTD.ŞTİBiological risk.