
REF	CONT	REF	CONT	REF	CONT
CHO-10100	2 x 50 mL	CHO-10300M	6 x 50 mL	CHO-10100A2	5 x 20 mL
CHO-10300	6 x 50 mL	CHO-10300P	6 x 50 mL		
CHO-10600	6 x 100 mL	CHO-10240M2	6 x 40 mL		
CHO-10180A	4 x 45 mL	CHO-10240M3	6 x 40 mL		

**KULLANIM AMACI**

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

INTENDED USE

For the *in vitro* quantitative determination of Total Cholesterol in serum

METHODOLOGY

Alain and Roeschlau began using cholesterol esterase and oxidase, in a single reagent to determine total cholesterol in serum. Trinder's color system of peroxidase / phenol/4-aminoantipyrine has been used successfully for some time now. The system's only drawback was the corrosive properties of phenol. The present method utilizes a phenol substitute that performs like phenol but without being corrosive.

Principle:

The intensity of the red color produced is directly proportional to the total cholesterol in the sample when read at 510(500±20)nm

REAGENT COMPOSITION

Active Ingredients	Concentration
PIPES	90 mM
4-Aminoantipyrin	0.40 mM
Cholesterol oxidase	350 U/L
Cholesterol esterase	500 U/L
Peroxidase	1000 U/L
Phenol	28 mM
pH	6,8 ±0.2

PRECAUTIONS

- This reagent is for *in vitro* diagnostic use only.
- Reagent contains sodium azide. Poison. Do not ingest.

REAGENT PREPARATION

Reagent is supplied ready to use.

REAGENT STORAGE

Store the reagent at 2-8°C (refrigerated).

The reagent is stable until the expiration date when stored at 2-8°C.

REAGENT DETERIORATION

The reagent should not be used if:

- The reagent is turbid.
- The working reagent does not meet stated performance parameters.

SPECIMEN COLLECTION AND STORAGE

- Nonhemolyzed serum is recommended.
- Cholesterol in serum is reported stable for seven days at room temperature (18-25°C) and six months when frozen and properly protected against evaporation.

INTERFERENCES

Bilirubin: No interference up to 10 mg/dL.

Hemoglobin: No interference up to 300 mg/dL.

Lipemia: No interference in the presence of triglycerides up to 1500 mg/dL.

Ascorbic acid: No interference up to 5 mg/dL.

ADDITIONAL EQUIPMENT REQUIRED BUT NOT PROVIDED:

- A clinical chemistry analyzer capable maintaining constant temperature (37°C), and measuring absorbance at 510(500±20)nm.
- Iron-free deionized water and related equipment, e.g.: pipettes
- Analyzer specific consumables, e.g.: sample and read cups
- Control and calibrator materials

PROCEDURE

Wavelength	: 510(500±20) nm
Working temperature	: 37°C
Optical path	: 1 cm
Assay type	: Endpoint
Direction	: Increasing

	Blank	Standard	Sample
Reagent	1000 µL	1000 µL	1000 µL
Distilled water	10 µL	--	--
Standard	--	10 µL	--
Sample	--	--	10µL

Mix and then incubate for 5 minutes at 37°C. Measure the absorbance of sample and standard against the reagent blank.

PROCEDURE NOTES

- The reagent and sample volumes may be altered proportionally to accommodate various instrument requirements.
- Grossly lipemic serums require a "sample blank". Add 0.01ml (10µl) of sample to 1.0mL saline, mix and read the absorbance against water. Subtract this value from the patient absorbance to obtain the corrected reading.

CALCULATION

(Abs = Absorbance)

Abs. of sample x conc. of std. (mg/dL) = Cholesterol (mg/dL)

Abs. of standard

LIMITATIONS

Samples with values exceeding 700 mg/dL should be diluted 1:1 with saline, re-run and result multiplied by two.

CALIBRATION

Aqueous standards can be used to calibrate the procedure or an appropriate serum calibrator.

QUALITY CONTROL

The integrity of the reaction should be monitored by use of a two level control with known cholesterol values.

EXPECTED VALUES

Recommended Range:

Desirable Cholesterol: <200 mg/dL

Borderline-High Cholesterol: 200 – 239 mg/dL

High Cholesterol: >240 mg/dL

PERFORMANCE

Sensitivity/ Limit of Detection(LOD):

The lower limit of detection is 0,87 mg/dL.

Linearity:

When the recommendation is designed and studied, the test is linear up to 700 mg/dL. Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

Precision :

	Within run (n=20)		Between run (n=20)	
Mean (mg/dL)	87,5	160,1	88,0	159,5
SD	1,27	3,14	2,11	3,24
CV(%)	1,45	1,96	2,40	2,03

Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x). The results obtained using 40 samples were the following:

Correlation Coefficient (r) =0,9957

Regression $y = 1,0279x - 0,8277$

The results of the performance characteristics depend on the analyzer used.

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BİLİMSEL TIBBİ ÜRÜNLER PAZ.SAN.VE TİC.LTD.ŞTİ