REF	CONT	REF	CONT
CKMB-30050	1x40mL+1x10mL	CKMB-30150M	4x30mL+2x15mL
CKMB-30150	3x40mL+1x30mL	CKMB-30150P	4x30mL+2x15mL
CKMB-30160A	4x32mL+4x8mL	CKMB-30100M2	2x40mL+2x10mL
CKMB-30050A2	4x10mL+1x10mL	CKMB-30200M3	4x40mL+2x20mL

⌘ : 2-8 °C

### KULLANIM AMACI

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⚠ : 2-8 °C


### INTENDED USE

Quantitative determination of creatin kinase MB (CK-MB) in human serum .

### PRINCIPLE OF THE METOD

The procedure involves measurement of CK activity in the presence of an antibody to CK-M monomer. This antibody completely inhibits the activity of CK-MM and half of the activity of CK-MB while not affecting the B subunit activity of CK-MB and CK-BB. Then it's used the CK method to quantitatively determine CK-B activity. The CK-MB activity is obtained by multiplying the CK-B activity by two.

### CLINICAL SIGNIFICANCE

CK-MB is an enzyme formed by the association of two subunits from muscle (M) and nerve cells (B). CK-MB is usually present in serum at low concentration; it is increased after an acute infarct of myocardium and later descends at normal levels. Also is increased, rarely, in skeletal muscle damage. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data

### REAGENT COMPOSITION

Active Ingredients	Concentration
<b>REAGENT 1</b>	
Imidazol buffer pH 6.7	125 mM
D-Glucose	25 mM
NAC	25 mM
Magnesium acetate	12,5 mM
NADP	2,52 mM
Hexokinase (microbial)	≥6800U/L U/L
EDTA	2,02 mM
Anti-human polyclonal CK-M antibody (sheep)inhibition capacity	2000 U/L
<b>REAGENT 2</b>	
ADP	15,2 mM
AMP	25 mM
Di(adenozin) pentaphosphate	103 mM
G6PDH (microbial)	≥8800U/L
Creatine phosphate	250 mM

### REAGENT PREPARATION

Reagents are supplied in a two vial, ready to use, liquid form. For some analyzers, the reagents can be combined to make a working solution by mixing 4 part of Reagent 1 with 1 part of Reagent 2 (e.g., 20mL R1 + 5mL R2).

### Precautions:

Reagent 1 and Reagent 2:  
H360d May damage fertility or the unborn child  
P308+313 If exposed or concerned: get medical advice /attention

### REAGENT STORAGE & STABILITY

- Store the reagent at 2-8°C (refrigerated).
- The reagent is stable until the expiration date when stored at 2-8°C.
- Working reagent is stable for 7days at 2-8°C or 12hours at 20-25°C. Protect from light.

### REAGENT DETERIORATION

Do not use the reagent if:  
the reagent blank absorbance < 0,300 at 340 nm against water.

### SPECIMEN COLLECTION AND STORAGE

- Serum free of hemolysis or heparin plasma.
- Avoid exposure of samples to strong light.
- Serum is stable for 7 days at (2-8°C)
- CK-MB activity decreases a 10% after 24 hours at 4°C or 1 hour at 25°C.

### INTERFERENCES

**Lipemia:** No interference in the presence of triglycerides up to 1500 mg/dL.  
**Ascorbic acid:** No interference up to 20 mg/dL.

### ADDITIONAL EQUIPMENT REQUIRED BUT NOT PROVIDED

- A clinical chemistry analyzer capable maintaining constant temperature 37°C, and measuring absorbance at 340nm.
- Deionized water and related equipment, e.g.: pipettes
- Analyzer specific consumables, e.g.: sample cups
- Control materials.

### PROCEDURE

Wavelength	: 340 nm
Temperature	: 37°C
Optical path	: 1cm
Assay type	: Kinetic
Direction	: Increasing

- Bring reagent to room temperature (15-30°C).

- Set the photometer to 0 (zero) absorbance with distilled water.

	Sample
Reagent 1	800 µL
Reagent 2	200 µL
Sample	40 µL

- Mix and incubate 10min at 37° C.
- Read the record initial absorbance(A1).
- Repeat the absorbance reading exactly after 5 minutes(A2).
- Calculate the average change in absorbance . (ΔAbs=A2-A1)

### CALCULATIONS

Abs=absorbans

ΔAbs x 1651 = CK-MB Activity (U/L)

ΔAbs x 825 = CK-B Activity (U/L)

Calculating factor in automatic analyzers by kinetic method( ΔAbs/min ) is 8255.

### Percentage of CK-MB Activity in Sample

This is calculated by dividing the found value of CK-MB activity by the total CK activity and multiplying the ratio by 100 as follows:

$\frac{\text{CK-MB Activity}}{\text{Total CK Activity}} \times 100 = \% \text{ CK-MB Activity}$

Total CK Activity

**Temperature conversion factors:**To correct results to other temperatures multiply by:

Assay temperature	Conversion factor to		
	25°C	30°C	37°C
25°C	1,00	1,53	2,38
30°C	0,65	1,00	1,56
37°C	0,42	0,64	1,00

### LIMITATIONS

- If the assay of the total CK as reported above is greater than 2000 U/L at 37°C dilute the sample appropriately with physiological saline before assay of CK-MB. Multiply the result by the dilution factor to obtain the correct value of the isoenzyme.
- The method will also measure any CK-BB isoenzyme present in serum. The activity of this isoenzyme is usually negligible, however, if a significant amount of CK-BB activity is present the CK-MB activity will be overestimated.
- A macro form of BB (immunoglobulin complexed) has been observed which will be measured as B in this assay. If the measured CK-B activity is greater than 20% of the total CK activity the presence of Macro BB should be suspected.

### QUALITY CONTROL

The integrity of the reaction should be monitored by use of a two level control with known CK-MB values.

### EXPECTED VALUES

	25°C	30°C	37°C
<b>CK-MB</b>	<b>&gt; 10 U/L</b>	<b>&gt; 15 U/L</b>	<b>&gt; 24 U/L</b>
<b>TOTAL CK</b>	<b>25°C</b>	<b>30°C</b>	<b>37°C</b>
erkek	80 U/L	130 U/L	195 U/L
kadin	70 U/L	110 U/L	170 U/L

### PERFORMANCE

#### Sensitivity/ Limit of Detection(LOD):

The lower limit of detection is 1,67 U/L.

#### Linearity:

When the recommendation is designed and studied, the test is linear up to 350 U/L. Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

#### Precision :

	Within run (n=20)	Between run (n=20)	
Mean (U/L)	39,4	37,3	92,1
SD	1,10	1,90	1,96
CV(%)	2,79	1,94	2,13

#### Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x).The results obtained using 40 samples were the following:

Correlation Coefficient (r) =0,9998

Regression y=1,001x +5,5164

The results of the performance characteristics depend on the analyzer used.

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BİLİMSEL TIBBİ ÜRÜNLER PAZ.SAN.VE TİC.LTD.ŞTİ