KULLANIM AMACI

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INTENDED USE

For the quantitative determination of copper in serum or plasma. For in vitro diagnostic use only.

METHODOLOGY

Serum copper is routinely determined by atomic absorption spectrophotometry (AAS) or spectrophotometry in clinical laboratories. AAS is the preferred method but these methods generally involve the reaction of a protein-free filtrate with organic solvents or reagents for deproteinizing, and none is satisfactory in terms of sensitivity, rapidity, and easy of measurement. In this method compounds are water soluble, and react easily and quickly with the metals forming intensely colored complexes.

Principle:

In an acid solution (pH 4.7), copper is released from The Ceruloplasmin protein and reduced. Afterwards the cuprous ion reacts with the 3,5-DiBr-PAESA.[4-(3,5-dibromo-2-pyridylazo)-N-ethyl-N-(3-sulfopropyl) aniline] stain and forms a stable coloured complex which is photometrically measured at 580 nm.

The colour intensity of this complex is proportional to the amount of copper in the sample.

REAGENT COMPOSITION

Acetate buffer pH 4,7	≤ 0,022 mM
3,5-DiBr-PAESA	0,2 mM
Surfactant	≤ % 1
Preservatives and stabilizers	

REAGENT PREPARATION

Reagent is ready to use liquid.

REAGENT STORAGE

Up to the expiration date at 2°C to 8°C
28 days of stability after opening the vial.

PRECAUTIONS

1. This reagent is for *in vitro* diagnostic use only.
2. All specimens used in the test should be considered potentially infectious.

SPECIMEN COLLECTION AND STORAGE

1. Serum or plasma. Use only heparin anticoagulants. Stability: 24 hours at 2-8°C or 15 days at -20°C.
2. When stored at 2-8°C components of this kit will remain stable until the expiration date stated on the label.

INTERFERENCES

No interferences were observed to bilirubin up to 15 mg/dL , hemoglobin up to 0,5 mg/dL and triglycerides up to 1000 mg/dL.

MATERIALS REQUIRED BUT NOT PROVIDED

1. Clinical chemistry analyzer or colorimeter capable maintaining constant temperature (37°C) and measuring absorbans at 580 nm
2. Deionize water and related equipment e.g. pipettes
3. Analyzer specific consumables, e.g: sample and read cups
4. Control and calibrator materials

PROCEDURE

Wavelength	: 580 nm
Working temperature	: 37°C
Optical path	: 1 cm
Assay type	: Endpoint
Direction	: Increasing

	Blank	Standard	Sample
Reagent	1000 µL	1000 µL	1000 µL
Standard	--	50 µL	--
Sample	--	--	50 µL

Incubate 5 minutes at 37°C after mix. Read sample and standard absorbans against blank.

CALCULATION

(Abs=absorbance)

$\text{Abs. Sample} \times \text{Concentration standard} = \text{COPPER } (\mu\text{g/dL})$

Abs. Std.

Conversion factor : $(\mu\text{g/dL}) \times 0,157 = \mu\text{mol/L}$

CALIBRATION

bt products CALB-30 chemistry calibrator or COPPER standard must be used.

QUALITY CONTROL

The validity of the reaction should be monitored by use of the control sera with known normal and abnormal copper values.

EXPECTED VALUES

MEN 70 – 140 µg/dL

WOMEN 80 – 155 µg/dL

These values should only be used as a guideline.

Each laboratory should establish its own reference range.

PERFORMANCE

1. **Linearity limit:** Up to 500 µg/dL(78,5 µmol/L). Samples with copper concentration higher than 500 µg/dL must be diluted 1:10 with % 0.9 NaCl , repeat test and multiply the result by 10.

2. **Detection Limit :** Values less than 3 µg/dL give non-reproducible results.

3. **Precision:** N= 20 replicates

	Within Run (n= 20)			Run to Run (n=20)		
Mean(µg/dL)	72.6	121.2	170.0	72.6	121.2	170.0
SD	2.41	2.98	1.97	0.96	1.06	1.76
CV(%)	3.05	2.19	1.57	3.32	2.46	1.16

4. **Accuracy:** Comparison between this method(y) and another commercial method (x) gave the following results:

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BİLİMSEL TIBBİ ÜRÜNLER PAZ.SAN.VE TİC.LTD.ŞTİ