
REF	CONT	REF	CONT
CRE-10040	1x30mL+1x10mL	CRE-10240P	4x45mL+2x30mL
CRE-10200	3x50mL+1x50mL	CRE-10160M2	4x30mL+2x20mL
CRE-10160A	4x30mL+4x10mL	CRE-10160M3	3x40mL+2x20mL
CRE-10240M	4x45mL+2x30mL	CRE-10100A2	5x15mL+2x12,5mL

KULLANIM AMACI

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
INTENDED USE

For the *in vitro* quantitative determination of Creatinine in serum and urine

METHODOLOGY

In the first reaction, creatinase and sarcosine oxidase were used in the enzymatic hydrolysis of endogenous creatine to produce hydrogen peroxide, that is eliminated by catalase. In the second reaction, the catalase is inhibited by sodium azide and creatinase and 4- aminoantipyrine (4-AAP) were added, and only the creatine generated from creatinine by creatininase was hydrolyzed sequentially by creatinase and sarcosine oxidase to produce hydrogen peroxide. This newly-formed hydrogen peroxide was measured in a coupled reaction catalyzed by peroxidase, with N-ethyl-n-sulphopropyl-m-toluidine (TOPS)/4-AAP as a chromogen.

REAGENT COMPOSITION

Active Ingredients

REAGENT 1	Concentration
MOPS	25 mM
TOPS	0.5 mM
Creatinase	10 KU/L
Sarcosine oxidase	5 KU/L
Catalase	3 KU/L
EDTA	1 Mm
pH 7,5 ± 0.2	
REAGENT 2	
MOPS	90 mM
Creatininase	30 KU/L
Peroxidase	
Sodium azide	0.95 g/L
pH 8,0 ± 0.2	

REAGENT PREPARATION

Reagents are supplied in a two vial, ready to use, liquid form.

PRECAUTIONS

Reagent 1 contains MOPS ; R36/37/38 Irritation to eyes, respiratory system and skin. Keep locked up, Keep out of the reach of children.

In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective cloth gloves and eye/face protection.

REAGENT STORAGE

1. Store the reagent at 2 - 8°C .
2. All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2 - 8°C , protected from light and contamination.
3. R1 and R2 are stable 8 weeks after opening bottle.

REAGENT DETERIORATION

The reagent should not be used if:

1. The reagent is cloudy (contaminated).
2. The reagent fails to meet stated parameters of performance.

SPECIMEN COLLECTION AND STORAGE

1. Unhemolyzed serum and plasma.
2. Urine should be diluted 1: 50 with distilled water.
3. Creatinine is stable for 1 day at (2-8°C) .

INTERFERENCE

Bilirubin: No interference up to 40 mg/dL.

Hemoglobin: No interference up to 500 mg/dL.

Lipemia: No interference in the presence of triglycerides up to 1500 mg/dL.

Ascorbic acid: No interference up to 25 mg/dL.

Some drugs and substances may interfere accuracy of the creatinine.

ADDITIONAL EQUIPMENT REQUIRED BUT NOT PROVIDED

1. A clinical chemistry analyzer capable maintaining constant temperature (37°C), and measuring absorbance at 545 nm (525 -565 nm).
2. Deionized water and related equipment, e.g.: pipettes
3. Analyzer specific consumables, e.g.: sample and read cups
4. Control and calibrator materials.

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Rev.Tarih / No: 30.08.2023 / 7

1. Bring the reagents and instrument to 37°C.

2. Pipette into cuvette:

	Blank	Standard	Sample
Reagent1 (R1)	225 µL	225 µL	225 µL
Distilled water	5 µL	--	--
Standard	--	5 µL	--
Sample	--	--	5 µL

3. Mix and incubate at 37°C, 5 minutes

4. Read the absorbance (Abs.1) of the standard or calibrator and the samples, at 545 nm against the blank.

5. Add R2:

	Blank	Standard	Sample
Reagent2 (R2)	75 µL	75 µL	75 µL

6. Mix and incubate at 37°C ,5 minutes.

7. Read the absorbance (Abs.2) of the standard or calibrator and the samples, at 545 nm against the blank

CALCULATIONS:

(Abs = Absorbance)

$$\frac{\Delta \text{Abs Sample} - \Delta \text{Abs Blank}}{\Delta \text{Abs standard} - \Delta \text{Abs Blank}} \times \text{Conc. of std (mg/dL)} = \text{Creatinine (mg/dL)}$$

LIMITATIONS

Samples with values exceeding 25 mg/dL should be diluted 1:2 with saline and re-run. The final answer should be multiplied by two.

CALIBRATION

Use an aqueous creatinine standard or an appropriate serum calibrator.

QUALITY CONTROL

The integrity of the reagent should be monitored by use of a two level control with known Creatinine values.

EXPECTED VALUE

Serum or plasma	:Men	0.9 – 1.3 mg/dL
	Women	0,6 – 1.1 mg/dL
Urine	: Men	14 – 26 mg/kg/24h
	Women	11 – 20 mg/kg/24h

It is highly recommended that each laboratory establish its own reference range.

PERFORMANCE

Sensitivity/ Limit of Detection (LOD):

The lower limit of detection is 0,086 mg/dL.

Linearity:

When the recommendation is designed and studied, the test is linear up to 25 mg/dL. Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

Precision :

	Within run (n=20)		Between run (n=20)	
Mean (mg/dL)	1,02	3,88	1,01	3,95
SD	0,05	0,07	0,02	0,05
CV(%)	4,75	1,74	2,05	1,32

Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x).The results obtained using 40 samples were the following:

Correlation Coefficient (r) =0,9968

Regression y =0,9985 x + 0,1102

The results of the performance characteristics depend on the analyzer used

REFERENCES

1. Fassati et al.Clin.Chem.1983;29:1494-1496
2. Tietz Textbook of Clinical Chemistry, 3 rd edt.Burtis CA,Ashwood ER.WB Saunders Co.1999.
3. Young, D.S. effects of drugs on clinical lab.test 4th edt.AACC Press.1995.
4. Young, D.S. effects of drugs on clinical lab.test 4th edt.AACC Press.2001