KULLANIM AMACI

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Rev.Date / No: 30.08.2023 / 12

INTENDED USE

 For the *in vitro* quantitative determination of Creatinine in serum and urine

METHODOLOGY

Jaffe described a method in 1886 for the determination of creatinine involving a protein free filtrate and a reaction with picric acid in alkaline solution. Although since then several methods have been described the classic Jaffe reaction method is still the most widely used. The Jaffe reaction is subject to interferences by a number of substances, including protein and glucose. Modifications of the procedure have been developed to combat the drawbacks. The kinetic procedures have become popular because they are fast, simple and avoid interferences. This method is based on a modification of the above procedure, incorporating a surfactant and other ingredients to minimize protein and carbohydrate interferences.

Principle: Alkali Medium

 Creatinine + Sodium Picrate $\xrightarrow{\text{Alkali Medium}}$ Creatinine-picrate complex (yellow-orange)

Creatinine reacts with picric acid in alkaline conditions to form a color complex which absorbs at 510(505-515) nm. The rate of formation of color is proportional to the creatinine in the sample.

REAGENT COMPOSITION
Active Ingredients
REAGENT 1

 Sodium hydroxide
 pH 13.0 ± 0.2

Concentration

350 mM

REAGENT 2

 Picric acid
 pH 2,6 ± 0.1

10 mM

Precautions

1. This reagent is for *in vitro* diagnostic use only.
2. Picric Acid is a strong oxidizing agent. Avoid contact with skin
3. Sodium hydroxide is an alkali. Avoid ingestion and contact.

Reagent 1 and Reagent 2:
H314 Causes severe skin burns and serious eye damage.
P303+P361+P353 IF ON SKIN: Remove/take off immediately all contaminated clothing. Rinse skin with water/shower
P304+P340 IF INHALED: Remove victim to fresh air and keep at rest in position comfortable for breathing.
P305+P351+P338 IF IN EYES: Rinse continuously with water for several minutes.
Remove contact lenses if present and easy to do. Continue rinsing. Wear suitable gloves and eye/face protection.
REAGENT PREPARATION

Reagents are supplied in a two vial, ready to use, liquid form.

For example analyzers, the reagents can be combined to make a working solution by mixing 1 part of Reagent 1 with 1 part of Reagent 2. (e.g., 5 mL R1 + 5 mL R2)

REAGENT STORAGE

1. Store the reagent at 15-25°C (room temperature).
2. The reagent is stable until the expiration date when stored at 15-25°C.
3. Working reagent is stable for 24 hours when stored at 15-25°C.

REAGENT DETERIORATION

The reagent should not be used if:

1. The reagent is cloudy (contaminated).
2. The reagent fails to meet stated parameters of performance.

SPECIMEN COLLECTION AND STORAGE

1. Unhemolyzed serum is recommended.
2. Urine should be diluted to a final concentration of approximately 34 to 64 mg/dL. 1:100 dilution is recommended.
3. Creatinine in serum is stable for 7 days at refrigerated temperatures (-2-8°C) and for several months when frozen (-20°C) and protected from evaporation and contamination. Urine samples are stable for at least 7 days at 4°C.

INTERFERENCES
Hemoglobin: No interference up to 500 mg/dL.

Lipemia: No interference in the presence of triglycerides up to 1500 mg/dL.

Ascorbic acid: No interference up to 25 mg/dL.

ADDITIONAL EQUIPMENT REQUIRED BUT NOT PROVIDED

1. A clinical chemistry analyzer capable maintaining constant temperature (37°C), and measuring absorbance at 510 (505-515nm).
2. Deionized water and related equipment, e.g.: pipettes
3. Analyzer specific consumables, e.g.: sample and read cups

4. Control and calibrator materials.
PROCEDURE

 Wavelength : 510 (505-515) nm
 Working temperature : 37°C
 Optical path : 1 cm
 Assay type : Fixed time
 Direction : Increase

1. Set to photometer to 0 (zero) absorbance with distilled water.

	Standard	Sample
Reagent 1 (R1)	500 µL	500 µL
Reagent 2 (R2)	500 µL	500 µL
Standard	100 µL	--
Sample	--	100 µL

2. Mix and incubate at 37°C.
3. After 60 sec., read the record absorbance. (Abs.1)
4. Return tube to 37°C and after 120 sec. read the final absorbance (Abs.2)
5. Determine the ΔAbs. for sample and standard.

Procedure Note:

The reagent and sample volumes may be altered proportionally to accommodate various instrument requirements.

CALCULATIONS:

(Abs = Absorbance)

 $\frac{\Delta \text{Abs Sample}}{\Delta \text{Abs standard}} \times \text{Concentration of standard (mg/dL)} = \text{Creatinine (mg/dL)}$
LIMITATIONS

Samples with values exceeding 25 mg/dL should be diluted 1:1 with saline and re-run. The final answer should be multiplied by two.

CALIBRATION

Use an aqueous creatinine standard or an appropriate serum calibrator.

QUALITY CONTROL

The integrity of the reagent should be monitored by use of a two level control with known Creatinine values.

EXPECTED VALUE

Serum: 0.60 – 1.40 mg/dl

Urine: 0.80 – 2.00 g/day

It is highly recommended that each laboratory establish its own reference range.

PERFORMANCE
Sensitivity/ Limit of Detection(LOD):

The lower limit of detection is 0,014 mg/dL.

Linearity:

When the recommendation is designed and studied, the test is linear up to 25 mg/dL. Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

Precision :

	Within run (n=20)		Between run (n=20)	
Mean (mg/dL)	1,13	3,52	1,11	3,21
SD	0,02	0,02	0,04	0,21
CV(%)	1,53	0,056	3,19	3,82

Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x). The results obtained using 40 samples were the following:

Correlation Coefficient (r) = 0,9998

 Regression $y = 1,0091x - 0,0708$

The results of the performance characteristics depend on the analyzer used.

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1. Tilzer L.L., Jacobs D.S. "Laboratory Test Handbook" 3rd Ed, Lexi-Comp Inc., p.202 (1994).
BİLİMSEL TIBBİ ÜRÜNLER PAZ.SAN.VE TİC.LTD.ŞTİ