### KULLANIM AMACI

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
### INTENDED USE

For quantitative determination of ferritin in human serum.

### Principle:

The latex particles coated with anti human ferritin are agglutinated when they react with samples that contain ferritin. The latex particles agglutination is proportional to the concentration of the ferritin in the sample and can be measured by turbidimetry.

### REAGENTS

**R1 (DILUENT):** Glycine buffer 20 mmol/L, pH 8.5

**R2 (LATEX):** Latex particles coted with polyclonal anti human ferritin. pH 8.2

### CALIBRATION (Optional)

Serum Ferritin Calibrator is recommended.

Ferritin concentration is stated on the label vial.

### PRECAUTIONS

1. The reagents contain sodium azide 0.95 g/L.. Avoid any contact with skin or mucous.
2. The reagents from human donors have been given negative results to anti HIV 1/2, HbsAg and anti HCV . Handle cautiously is recommended.

### REAGENT PREPARATION

R1 :Ready to use.

R2 :Ready to use. Shake gently the vial before use.

### STORAGE AND STABILITY

1. The reagents will remain stable until the expiration date printed on the label, when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use the reagents after the expiration date.
2. Protect reagents from extreme heat, light or freezing.

### REAGENT DETERIORATION

1. Increment of blank reagent absorbance.
2. Presence of particles, turbidity.

### SAMPLES

1. Fresh serum. Stable 7 days at 2-8°C, 3 months at - 20°C.
2. Samples with presence of fibrin should be centrifuged before testing.
3. Hemolyzed or contaminated samples are not suitable for testing.

### INTERFERENCES

Bilirubin (20 mg/dL), hemoglobin (10 g/L), rheumatoid factor (600 IU/mL), do not interfere. Lipemia interferes.

### ADDITIONAL EQUIPMENT

Spectrophotometer or photometer thermostable at 37°C with a (540±20) nm filter.

### PROCEDURE

Wavelength : 540 nm (530-550 nm)  
 Working temperature : 37°C  
 Optical path : 1 cm  
 Assay type : Turbidimetry  
 Direction : Increasing

1. Bring the reagents and the photometer (cuvette holder) to 37°C .
2. Adjust the instrument to zero with distilled water.
3. Pipette into a test tube

Sample or calibrator tube	
R1 (DILUENT)	400 µL
Calibrator or sample	45 µL

4. Mix and pipette R2 Reagent into cuvette

R2 (LATEX)	100 µL
------------	--------

5. Mix well and record the absorbance immediately (Abs.1) and after 5 minutes of the reagent R2 addition read the final absorbance (Abs. 2)

### CALCULATION

Calculate the absorbance difference (Abs.2-Abs.1) of each point of the calibration curve and plot the values obtained against the ferritin concentration of each calibrator dilution. Ferritin concentration in the sample is calculated by interpolation of its(Abs.2-Abs.1) in the calibration curve.

### QUALITY CONTROL

Ferritin control sera are recommended to monitor the performance of manual and automated assay procedures. Each laboratory should establish its own Quality control scheme and corrective actions if controls do not meet the acceptable tolerances.

### EXPECTED VALUES

Children : 7 – 140 µg/L

Men : 20 – 250 µg/L

Women : 20 – 200 µg/L

Each laboratory should establish its own reference range.

### CALIBRATION

Prepare dilutions of the calibrator using NaCl 0.9 g/dL as diluent. Multiply the concentration of the calibrator by the corresponding factor indicated in the table below to obtain the ferritin concentration of each point of the curve.

DILUTION	1	2	3	4	5
Ferritin Calibrator (µL)	--	25	50	75	100
NaCl 0.9 g/dL( µL)	100	75	50	25	--
Factor	0.0	0.25	0.5	0.75	1.0

### PERFORMANCE

**1.Linearity:** Up to 300 µg/L, under the described assay conditions.Samples with higher concentrations should be diluted 1/5 with NaCl 0.9 g/dL and retest again.

**2.Detection Limit:**

Values less than 3 µg/L give non-reproducible results.

**3.Sensitivity:** 2.7 mA/ µg/L.

**4.Prozone effect :** Up to 4000 µg/L.

**5.Precision:**

	intra assay(n=10)		inter-assay(n=10)	
Mean(µg/L)	65	178	65	178
CV%	3.56	1.87	5.16	2.9

**6. Accuracy:** Results obtained using this reagent did not show differences when compared with other commercial reagents . Studies of comparison are available upon request.

### NOTES :

1. Calibrator dilutions in plastic tubes should be avoided. Ferritin antigen may coat to the walls.
2. This method may be used with different instruments .Any application to an instrument should be validated to demonstrate that results meets the performance characteristics of the method. It is recommended to validate periodically the instrument. Contact to the distributor of any question on the application method.
3. The linearity limit depends on the sample/reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
4. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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Bilimsel Tıbbi Ürünler Paz.San.Tic.Ltd.Şti.