
REF	CONT	REF	CONT	REF	CONT
GGT-10100	2x40mL+2x10mL	GGT-10228M	4x45mL+2x24mL	GGT-10100A2	4x20mL+2x10mL
GGT-10250	5x40mL+1x50mL	GGT-10228P	4x45mL+2x24mL		
GGT-10500	5x80mL+1x100mL	GGT-10200M2	4x40mL+2x20mL		
GGT-10160A	4x32mL+4x8mL	GGT-10200M3	4x40mL+2x20mL		

KULLANIM AMACI

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BİLİMSEL KAYNAKÇA

- Gendler S.-GT, Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1120-1123.
- Persijn J P et al. J Clin Chem Clin Biochem 1976; (14) 9: 421-427.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACCC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACCC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACCC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACCC 1995.

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
INTENDED USE

This reagent is intended for the *in vitro* quantitative kinetic determination of gamma glutamyltransferase (Gamma-GT) in human serum.

METHODOLOGY

Methods for determining Gamma-GT are based on the use of glutamyl derivatives of aromatic amines as substrate material. Gamma-GT reagent uses this soluble 3-carboxyl derivative. This method is based on the kinetic photometric test, according to the International Federation of Clinical Chemistry and Laboratory Medicine (IFCC).

Principle:
 L-Gamma-Glutamyl-3-carboxy-4-nitroanilide + Glycylglycine \longrightarrow L-Gamma-glutamyl-glycylglycine + 2-nitro 5 aminobenzoic acid

Gamma-GT in the sample catalyzes the transfer of the glutamyl group from L-Gamma-glutamyl-3-carboxy-4-nitroanilide to glycylglycine according to the above reaction. The amount of 2-nitro 5 aminobenzoic acid formed is proportional to Gamma-GT activity and may be measured kinetically at 405 nm by the increasing intensity of the yellow color formed.

REAGENT COMPOSITION

Active Ingredients

REAGENT 1	Concentration
Tris buffer	100 mM
Glycylglycine	100 mM
pH 8,55± 0,1	
REAKTİF 2	Concentration
L-gamma glutamyl 3-carboxy 4-nitroanilide	3.0 mM
pH 6,15± 0,1	

Precautions and Warnings:

- For *in vitro* diagnostic use only.
- DO NOT pipette by mouth. Avoid contact with skin and eyes. If spilt, thoroughly wash affected area with water
- Reagent contains Sodium Azide as a preservative. Do not use the reagent after the expiration date printed on the kit.

REAGENT PREPARATION

Reagents are ready to use, liquid form.
 Working Reagent : The reagents can be combined to make a working solution by mixing 4 parts of Reagent 1 with 1 part of Reagent 2 (e.g., 20mL Rgt 1 to 5 mL Rgt 2).

STABILITY AND STORAGE

- Store the reagents at 2-8°C (refrigerated).
- The reagents are stable until the expiration date when stored at 2-8°C.
- Working reagent is stable for 21 days when stored tightly closed at (2-8°C), or 5 days at room temperature protected from light and contamination during their use.
- Do not freeze the reagents.
- Reagent 2 should be protected from light.

REAGENT DETERIORATION

DO NOT USE THE REAGENT IF:
 1. The reagent is turbid and presence particles.
 2. The reagent has an optical density greater than 1.2 at 405 nm.

SPECIMEN COLLECTION AND STORAGE

- Use non-hemolyzed serum
- Serum samples should be stored 3 days at 2-8°C, 8 hours at(15-25°C) and 1 month at (-20°C).

INTERFERENCES

Bilirubin: No interference up to 40 mg/dL.
Hemoglobin: No interference up to 100 mg/dL.
Lipemia: No interference in the presence of triglycerides up to 1500 mg/dL.
Ascorbic acid: No interference up to 25 mg/dL

ADDITIONAL EQUIPMENT REQUIRED BUT NOT PROVIDED

- A clinical chemistry analyzer capable maintaining constant temperature 25°C, 30°C or 37°C , and measuring absorbance at 405nm.
- Deionized water and related equipment, e.g.: pipettes
- Analyzer specific consumables, e.g.: read and sample cups
- Control materials.

PROCEDURE

Wavelength	: 405 nm
Working temperature	: 25°C /30°C /37°C
Optical path	: 1cm
Assay type	: Kinetic
Direction	: Increasing

- Bring reagent to room temperature (15-30°C).
- Set the photometer to 0 (zero) absorbance with distilled water.

Reagent 1	Sample
Reagent 2	800 µL
Sample	200 µL
	100 µL

- Mix and incubate at 37° C.
- After 60 sec., read and record absorbance.
- Return tube to 37° C.
- Repeat readings every minute for the next three minutes
- Calculate the average absorbance difference per minute. (ΔAbs/Min.)

CALCULATIONS

One Unit (U/L) is defined as the amount of enzyme that catalyzes the transformation of one micromole of substrate per minute under the specified conditions.
 GGT U/L = Δ Abs./min. x 1190

Temperature conversion factors

To correct results to other temperatures multiply by:

Assay temperatures	25°C	30°C	37°C
25°C	1.00	1.37	1.79
30°C	0.73	1.00	1.30
37°C	0.56	0.77	1.00

QUALITY CONTROL:

The integrity of the reagent should be monitored by use of a two level control with known GGT values.

EXPECTED VALUES

	25°C	30°C	37°C
Women	4-18 U/L	5-25 U/L	7-32 U/L
Men	6-28 U/L	8-38 U/L	11-50 U/L

It is strongly suggested that each laboratory establish its own normal range.

PERFORMANCE

Sensitivity/ Limit of Detection(LOD):

The lower limit of detection is 1,19 U/L.

Linearity:

When the recommendation is designed and studied, the test is linear up to 1200 U/L. Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

Precision :

	Within run (n=20)	Between run (n=20)
Mean (U/L)	47,6	218,2
SD	0,31	1,01
CV(%)	0,64	0,46
		0,51
		0,40

Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x).The results obtained using 40 samples were the following:
 Correlation Coefficient (r) =1,000

Regression $y=1,1524x +1,5314$

The results of the performance characteristics depend on the analyzer used.

REFERENCES

- Gendler S.-GT. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1120-1123.
- Persijn J P et al. J Clin Chem Clin Biochem 1976; (14) 9: 421-427.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AAC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AAC Press, 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AAC Press, 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AAC Press, 1995.

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BİLİMSEL TIBBİ ÜRÜNLER PAZ.SAN.VE TİC.LTD.ŞTİ