
REF	CONT	REF	CONT	REF	CONT
GLU -10100	2 x 50 mL	GLU-10300M	6 x 50 mL	GLU-10100A2	5 x 20 mL
GLU -10300	6 x 50 mL	GLU-10300P	6 x 50 mL		
GLU -10600	6 x 100 mL	GLU-10240M2	6 x 40 mL		
GLU-10180A	4 x 45 mL	GLU-10240M3	6 x 40 mL		

**KULLANIM AMACI**

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

INTENDED USE

For the in vitro quantitative determination of Glucose in serum.

METHODOLOGY

Early enzymatic methods for glucose determination used Glucose Oxidase to catalyze the oxidation of glucose to hydrogen peroxide and gluconic acid. The hydrogen peroxide that is formed is measured by the oxidation of a chromagen. Many chromagens were investigated but many were discarded because of possible carcinogenicity, toxicity, instability or because they were affected by many interfering substances. Trinder modified Emerson to develop an efficient peroxidase phenolaminophenazone system for the quantitation of hydrogen peroxide by formulation of a red quinoneimine dye. This method is less influenced by interfering substances and does not suffer from the many drawbacks of earlier methods. The present procedure is based on the above principle but utilizes a non-corrosive phenol substitute for added safety and convenience.

PRINCIPLE:

REAGENT COMPOSITION

Glucose Oxidase	≤ 23 U/mL
Peroxidase	≤ 0,75 U/mL
Aminoantipyrene	0,30 mM
4-Chlorophenol	< 10 mM
Non-reactive stabilizers and fillers Sodium Azide 0,05%.	
pH 7.4 ± 0.15	

PRECAUTIONS

1. This reagent is for in vitro diagnostic use only.
2. This reagent contains sodium azide at 0,05%

REAGENT PREPARATION

Reagent is liquid and ready to use.

REAGENT STORAGE

1. The reagent should be stored refrigerated at 2-8°C.
2. The reagent is stable until the expiration date when stored at 2-8°C.
3. Do not freeze the reagent.

REAGENT DETERIORATION

Do not use if:

1. The reagent fails to recover stated control values or meet stated linearity.
2. The reagent develops turbidity, or other evidence of microbial growth.

SPECIMEN COLLECTION AND STORAGE

1. Non-hemolyzed serum or heparinized plasma is recommended.
2. Serum must be separated from the clot promptly since the rate of glucose decrease is approximately 7% per hour in whole blood
3. Glucose in serum or plasma is stable for 24 hours when stored refrigerated (2-8°C)

INTERFERENCES

Bilirubin: No interference up to 10 mg/dL.

Hemoglobin: No interference up to 500 mg/dL.

Lipemia: No interference in the presence of triglycerides up to 1500 mg/dL.

Ascorbic acid: No interference up to 10 mg/dL.

Grossly lipemic or icteric samples will cause false glucose values, consequently a patient blank should be run.

For a comprehensive list of interfering substances see Young, et al

ADDITIONAL EQUIPMENT REQUIRED BUT NOT PROVIDED

1. A clinical chemistry analyzer capable maintaining constant temperature (37°C), and measuring absorbance at 510nm.
2. Deionized water and related equipment, e.g.: pipettes.
3. Analyzer specific consumables, e.g.: sample and read cups.
4. Control and calibrator materials.

PROCEDURE

Wavelength	: 510 nm (500-550nm)
Working temperature	: 37°C
Optical path	: 1 cm
Assay type	: Endpoint
Direction	: Increasing

	Blank	Standard	Sample
Reagent	1000 µL	1000 µL	1000 µL
Distilled water	10 µL	--	--
Standard	--	10 µL	--
Sample	--	--	10 µL

Mix and then incubate for 10 minutes at 37°C. Measure the absorbance of sample and standart against the reagent blank.

Note:The reagent and sample volumes may be altered proportionally to accommodate various instrument requirements.

CALCULATIONS

(Abs = Absorbance)

$\frac{\text{Abs (sample)}}{\text{Abs (standard)}} \times \text{Con. of standard (mg/dL)} = \text{Glucose (mg/dL)}$

Abs (standard)

CONVERSION FACTOR

$\text{Glucose (mg/dL)} \times 0,05551 = \text{Glucose (mmol/L)}$

CALIBRATION

Use an aqueous Glucose standard or an appropriate serum calibrator.

QUALITY CONTROL

The integrity of the reaction should be monitored by use of a two level control with known Glucose values

EXPECTED VALUES

Normal range is reported to be: 70 -115 mg/dL.

It is strongly recommended that each laboratory establish its own normal range.

PERFORMANCE

Sensitivity/ Limit of Detection(LOD):

The lower limit of detection is 1,34 mg/dL.

Linearity:

When the recommendation is designed and studied, the test is linear up to 500 mg/dL. Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

Precision :

	Within run (n=20)		Between run (n=20)	
Mean (U/L)	106,2	233,9	106,2	235,7
SD	1,51	3,82	1,82	3,21
CV(%)	1,42	1,63	1,72	1,36

Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x).The results obtained using 40 samples were the following:

Correlation Coefficient (r) =0,9963

Regression $y=0,9228x +2,5493$

The results of the performance characteristics depend on the analyzer used.

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BİLİMSEL TIBBİ ÜRÜNLER PAZ.SAN.VE TİC.LTD.ŞTİ