REF	CONT
GLY-10040	1x30mL+1x10mL+Sam.Dil:2x50mL
GLY-10024A	1x18mL+1x6mL+Sam.Dil:2x30mL
GLY-10040A2	1x30mL+1x10mL+Sam.Dil:1x100mL
GLY-10080M	4x15mL+2x10mL+Sam.Dil:2x100mL

REF	CONT
GLY-10080P	4x15mL+2x10mL+Sam.Dil:2x100mL
GLY-10040M2	1x30mL+1x10mL+Sam.Dil:1x100mL
GLY-10080M3	2x30mL+2x10mL+Sam.Dil:2x100mL

: 2-8 °C

Numune, kalibratör,  
kontrol tüpü

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
**INTENDED USE**

Quantitative determination of Glycated hemoglobin(HbA1c) in human whole blood.

**METHODOLOGY**

This kit is based on latex agglutination. HbA1c in test samples is absorbed onto the surface of latex particles, which react with Anti-HbA1c (antigen-antibody reaction). The turbidity caused by latex agglutination is measured at 660 nm, and the HbA1c concentration in whole blood is calculated from a standard curve.

**REAGENTS**

R1 : Latex

R2 : Anti-HbA1c (Anti-human hemoglobin A1c mouse monoclonal antibody)

Sample Diluent : Distilled wWater and stabilizers.

**CALIBRATION (Optional)**

For only use **bt products** (GLYCAL-10) Auto-HbA1c Calibrator Five (5) Level.

**PRECAUTIONS**

1. For only in vitro diagnostic use.
2. Reagent R1 and R2 are light-sensitive. Store in a dark place.
3. Avoid ingestion and contact with skin and eyes.

**REAGENT PREPARATION**

All the reagents are ready to use.

**STORAGE AND STABILITY**

1. Reagents are stable until their expiration date when stored at 2-8°C.
2. Protect reagents from extreme heat, light or freezing.
3. One month at 2 - 10°C after opening

**REAGENT DETERIORATION**

1. Do not use the reagents after the expiration date.
2. Increment of blank reagent absorbance.
3. Presence of particles, turbidity.

**SAMPLE**

Venous whole blood samples collected with K2-EDTA anticoagulant can be used.

It is recommended that samples be used within 2 weeks of collection when stored refrigerated. Prior to testing, whole blood samples should be mixed by gentle inversion to re-suspend settled erythrocytes.

**INTERFERENCES**

High lipid (3000 formazin turbidity units) or bilirubin (conjugated type 30mg/dL, free type 30 mg/dL) counts in blood samples do not affect the test results.

**ADDITIONAL EQUIPMENT**

1. Spectrophotometer or photometer thermo stable at 37°C with a 660 nm filter.
2. (GLYCAL-10) Auto-HbA1c Calibrator Five (5) Level and (GLYCON-10) Auto- HbA1c Control are recommended.

**SAMPLE PRETREATMENT (HEMOLYSED SAMPLE)**

1. Dispense 500 µL of Sample Diluent in a sample cup.
2. Whole blood samples should be mixed by gentle inversion at least 5 times to resuspend settled erythrocytes.
3. Add 10 µL of fully resus-pended whole blood sample to the Sample Diluent in the sample cup.
4. Incubate at room temperature (25°C) for 5 min to completely lyse the red blood cells.

**PROCEDURE**

Wavelength: : 660 nm  
Working temperature: : 37°C  
Optical path : 1 cm  
Assay type: : Endpoint  
Direction: : Increasing

Sample or calibrator tube

Reconstituted R1	180 µL
Hemolysed Sample, Calibrator or Control	4,8 µL

1. Mix and incubate for 5 minutes at 37° C and read the absorbance immediately (Abs.1)

R2	60 µL
----	-------

2. Mix and incubate for 5 minutes at 37° C and read the absorbance immediately (Abs.2)

**CALIBRATION PROSEDURE**

Auto- HbA1c calibrator sets are intended for use with Auto- HbA1c assay reagents.

Calculate the absorbance difference (Abs.2-Abs.1) of each point of the calibration curve and plot the values obtained against the HbA1c concentration of each calibrator dilution. HbA1c concentration in the sample is calculated by interpolation of its(Abs.2-Abs.1) in the calibration curve.

**CALCULATIONS**

The HbA1c concentration is expressed directly as %HbA1c by use of a suitable calibration curve in which the calibrators have values for each level in %HbA1c.

No calculation step is needed.

**QUALITY CONTROL**

**bt products** GLYCON-10 Auto-HbA1c Controls are recommended to monitor the performance of manual and automated assay procedures. Each laboratory should establish its own Quality control scheme.

**EXPECTED VALUES**

NGSP (%) : 4,6 – 6,2 %  
IFCC (SI) mmol/mol) : 26,8 – 44,3 mmol/molHb

NGSP = (0,0915 X IFCC) + 2,152

Each laboratory should establish its own reference range to reflect the age,sex,diet and geographical location of the population.

**PERFORMANCE**

**1.Linearity** : The linearity of the assay is up to 13,5 0 % HbA1c. Samples with values above 13,5% should not be diluted and retested. Instead the values should be reported as higher than 13,5%

**2.Assay Range** : 3,3% - 13,5%.

**3.Correlation**

Study results in comparison with another latex agglutination method (inhouse

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Bilimsel Tıbbi Ürünler Paz.San.Tic.Ltd.Şti.