REF	CONT	REF	CONT	REF	CONT
GPT-11100	2x40mL+2x10mL	GPT-11228M	4x45mL+2x24mL	GPT-11100A2	4x20mL+2x10mL
GPT-11250	5x40mL+1x50mL	GPT-11228P	4x45mL+2x24mL		
GPT-11500	5x80mL+1x100mL	GPT-11300M2	6x40mL+3x20mL		
GPT-11160A	4x32mL+4x8mL	GPT-11200M3	4x40mL+2x24mL		

⚠ : 2-8 °C

KULLANIM AMACI

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⚠ : 2-8 °C



INTENDED USE

For the quantitative determination of Alanine Aminotransferase in serum.

METHODOLOGY

The International Federation of Clinical Chemistry (IFCC) published a proposed recommended method in 1980 utilizing the LDH-NADH coupled assay. The procedure described herein is based on that method.

Principle

ALT catalyzes the transfer of the amino group from L-alanine to α -ketoglutarate resulting in the formation of pyruvate and L-glutamate. Lactate dehydrogenase catalyzes the reduction of pyruvate and the simultaneous oxidation of NADH to NAD. The resulting rate of decrease in absorbance is directly proportional to ALT activity.

REAGENT COMPOSITION

Active Ingredients

Reagent 1	Concentration
Tris Buffer	125 mM
L-Alanine	680 mM
LDH (microbial)	>2000 U/L
pH 7.5 ± 0.1	
Reagent 2	Concentration
α -ketoglutaric acid	97 mM
NADH	1.1 mM
Sodium azide	0.01 %
pH 10.5 ± 0.1	

REAGENT PREPARATION

Prepare working reagent by mixing 4 parts of R1 reagent with 1 part R2 reagent. (e.g. 200 μ L R1 with 50 μ L R2 reagent.)

REAGENT STORAGE

- Store reagents at 2-8°C.
- Working reagent is stable for 5 days at room temp. (20-25°C) and for 2 weeks when refrigerated (2-8°C).

REAGENT DETERIORATION

Do not use reagent if:

- The initial absorbance at 340nm is below 0,8.
- The reagent fails to meet stated parameters of performance.

PRECAUTIONS

- This reagent set is for *in vitro* diagnostic use only.
- The R1 reagent contains sodium azide as a preservative. Do not ingest.

SPECIMEN COLLECTION AND STORAGE

- Non-hemolyzed serum is recommended. Red cells contain ALT which can give falsely elevated results.
- ALT in serum is stable for three days at room temperature (15-30°C), 1 week at 4°C.

INTERFERENCES

Bilirubin: No interference up to 40 mg/dL.
 Hemoglobin: No interference up to 500 mg/dL.
 Lipemia: No interference in the presence of triglycerides up to 1500 mg/dL.
 Ascorbic acid: No interference up to 25 mg/dL.

MATERIALS REQUIRED BUT NOT PROVIDED

- Accurate pipetting devices.
- Test tubes/rack.
- Timer.
- Spectrophotometer able to read at 340 nm. (UV)
- Heating bath/block (37°C).

PROCEDURE

Wavelength : 340 nm
 Assay Type : Kinetic
 Reaction Direction : Decreasing
 Temperature : 37°C
 Optical path : 1 cm

- Bring to room temperature (15 -30 °C)
- Set the photometer to 0 (zero) with distilled water

	Sample
Reagent 1	800 μ L
Reagent 2	200 μ L
Sample	100 μ L

- Mix and incubate 37 °C
- After 60 sec. read and record the absorbance
- Return tube to 37 °C.
- Repeat the absorbance reading exactly after 1,2 and 3 minutes.
- Calculate the average absorbance difference per minute. (Δ Abs/Min.) The Δ Abs/Min multiplied by the factor 1768 will yield results in U/L

Procedure Notes

- If the spectrophotometer being used is equipped with a temperature-controlled cuvette, the reaction mixture may be left in the cuvette while the absorbance readings are taken.
- A very low final reading, together with a small absorbance change between readings could indicate a very high ALT level. Dilute and re-assay as necessary.

LIMITATIONS

Turbid or highly icteric samples may give readings whose initial absorbance exceeds the capabilities of the spectrophotometer. More accurate results may be obtained by using 50 μ l of sample and multiplying the final answer by two.

CALIBRATION

The procedure is standardized by means of the millimolar absorptivity of NADH taken as 6.22 at 340nm under the test conditions described.

CALCULATION

One international Unit (U/L) is defined as the amount of enzyme that catalyzes the transformation of one micromole of substrate per minute under specified conditions.

$$\text{ALT (U/L)} = \Delta\text{Abs./min.} \times 1768$$

NOTE:

- If test parameters are altered the factor has to be recalculated using the above formula.
- SI Units: To convert to SI Units (nkat/L) multiply U/L by 16.67.

QUALITY CONTROL

The validity of the reaction should be monitored using control sera with known normal and abnormal ALT (SGPT) values.

EXPECTED VALUES

< 40 U/L at 37°C

It is strongly recommended that each laboratory establish its own normal range.

PERFORMANCE

Sensitivity/ Limit of Detection (LOD):

The lower limit of detection is 1,46 U/L.

Linearity :

When the recommendation is designed and studied, the test is linear up to 500 U/L . Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

Precision:

	Within run (n=20)		Between run (n=20)	
Mean (U/L)	45,6	119,7	44,5	116,4
SD	0,96	1,04	0,48	0,6
CV(%)	2,11	0,87	1,07	0,51

Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x). The results obtained using 40 samples were the following:

Correlation Coefficient (r) =0,9981

Regression equation $y=0,9983x -0,3206$ (Two reagents procedure)

The results of the performance characteristics depend on the analyzer used.

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BİLİMSEL TIBBİ ÜRÜNLER PAZ.SAN.VE TİC.LTD.ŞTİ