
REF	CONT	REF	CONT
IGG-30100	2x40mL+2x10mL	IGG-30050MP	2x20mL+1x10mL
IGG-30050A	4x10mL+2x5mL	IGG-30050M2	2x20mL+2x5mL
IGG-30040A2	4x8mL+1x8mL	IGG-30050M3	1x40mL+1x10mL

**INTENDED USE**

For quantitative determination of immunglobulin IgG in human serum or plasma.

**PRINCIPLE OF THE METHOD**

Anti human IgG antibodies when mixed with samples containing IgG form insoluble complexes. These complexes cause an absorbance change, dependent upon the IgG concentration of the patient sample, that can be quantified by comparison from a calibrator of known IgG concentration.

**REAGENTS**

**R1 (DILUENT):** Tris buffer 20 mmol/L, PEG 8000 , Ph 8.3.  
Sodium azide 0.95 g/L.

**R2 (ANTIBODY):** Goat serum, anti human IgG .Ph 7.5. Sodium azide 0.95 g/L

**CALIBRATION (Optional)**

Serum Proteins Calibrator is recommended.

**REAGENT PREPARATION**

Reagents are ready to use.

**STORAGE AND STABILITY**

- All the components of the kit are stable until the expiration date printed on the label, when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use the reagents after the expiration date.
- Do not freeze the reagents: frozen antibody or diluent could change the functionality of the test.

**REAGENT DETERIORATION**

Presence of particles and turbidity.

**SAMPLES**

- Fresh serum or plasma; EDTA or heparin should be used as anticoagulant.  
Stable 7 days at 2-8°C, 3 months at -20°C .
- Samples with presence of fibrin should be centrifuged before testing.
- Hemolyzed or lipemic samples are not suitable for testing.

**INTERFERENCES**

Bilirubin (20 mg/dL), hemoglobin (10.0 g/dL), lipemia (10 g/L) do not interfere. Rheumatoid factor may interfere at 300 IU/mL.

**ADDITIONAL EQUIPMENT**

- Spectrophotometer or photometer thermostable at 37°C with a 600 nm (580-620) nm filter.
- Serum Proteins Calibrator and Serum Proteins Control are recommended for manual and automated assay procedures.

**PROCEDURE**

Wavelength : 600 nm (580-620)  
Working temperature : 37°C  
Optical path : 1 cm  
Assay type : Turbidimetry  
Direction : Increasing

- Bring the reagents and the photometer (cuvette holder) to 37°C.
- Adjust the instrument to zero with distilled water.
- Pipette into a test tube.

	Sample or calibrator tube
<b>R1 (DILUENT)</b>	300 µL
Sample or calibrator	4 µL

- Mix and read the absorbance (Abs.1) after the sample or calibrator addition.
- Pipette into a test tube

<b>R2 (ANTIBODY)</b>	75 µL
----------------------	-------

- Mix well and read the absorbance (Abs.2) of calibrators and sample exactly 2 minutes after the R2 addition.

**CALCULATION**

Calculate the absorbance difference (A2-A1) of each point of the calibration curve and plot the values obtained against the IgG concentration of each calibrator dilution. IgG concentration in the sample is calculated by interpolation of its(A2-A1) in the calibration curve.

**QUALITY CONTROL**

Serum Proteins Control are recommended to monitor the performance of manual and automated assay procedures. Each laboratory should establish its own Quality control scheme and corrective actions if controls do not meet the acceptable tolerances.

**EXPECTED VALUES**

700 -1600 mg/dL.

Each laboratory should establish its own reference range.

**CALIBRATION**

Prepare dilutions of the calibrator using NaCl 0.9 g/dL as diluent. Multiply the concentration of the calibrator by the corresponding factor indicated in the table below to obtain the IgG concentration of each point of the curve.

DILUTION	1	2	3	4	5	6
CALIBRATORS (µL)	--	10	25	50	75	100
NaCl 0.9 g/dL(µL)	100	90	75	50	25	--
Factor	0.0	0.1	0.25	0.5	0.75	1.0

**PERFORMANCE**

- Linearity limit :** Up to 3000mg/dL, under the described assay conditions. Samples with higher concentrations should be diluted 1/5 with NaCl 0.9 g/dL and retest again. The linearity limit depends on the sample/reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Detection Limit:**  
Values less than 10.3 mg/dL give non-reproducible results.
- Analytical sensitivity:** Δ 0.6 mA / mg/dL at 359 mg/dL.
- Prozone effect:** Up to 8000 mg/dL.
- Precision:** The reagent has been tested for 20 days, using three levels of serum.

	CV%		
	340.3 mg/dL	801.96 mg/dL	1517.5 mg/dL
<b>Total</b>	2.1	2.8	4.8
<b>Within run</b>	0.9	0.7	1.0
<b>Between run</b>	1.5	1.5	1.8
<b>Between day</b>	1.0	2.2	4.4

- Accuracy:** Results obtained using this reagent(y) were compared to those obtained using another commercially available method (x) with similar characteristics. 79 samples ranging from 450 - 2600 mg/dl of IgG were assayed. The correlation coefficient( r ) was 0.94 and the regression equation  $y = 0.957x + 105.67$  .The result of the performance characteristics depend on the analyzer used.

**NOTES :**

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

**REFERENCES**

- Clinical guide to laboratory tests, edited by NW Tietz W B Saunders Co. Philadelphia 483, 1983
- Dati F et al. Eur. J Clin. Chem. Clin. Biochem 1996;34:517 - 520
- Pesce A J and Kpalian, LA, Methods in Clin. Chem. The CV Mosby Company, St. Louis MO 1987
- Kreutzer HJH, J Clin Chem Clin. Biochem 1976;14:401 -406
- Young DS. Effects of drugs on clinical laboratory tests. 4th ed. AACCC Pres. 1995
- Friedman and Young. Effects of disease on clinical lab. tests 3 rd ed. AACCC Pres. 1997

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**BİLİMSSEL KAYNAKÇA**

1. Clinical guide to laboratory tests, edited by NW Tietz W B Saunders Co. Philadelphia 483, 1983
2. Dati F et al. Eur. J Clin. Chem. Clin. Biochem 1996;34:517 - 520
3. Pesce A J and Kpalan, LA, Methods in Clin. Chem. The CV Mosby Company, St. Louis MO 1987
4. Kreutzer HJH, J Clin Chem Clin. Biochem 1976;14:401 -406
5. Young DS. Effects of drugs on clinical laboratory tests. 4th ed. AACC Pres. 1995
6. Friedman and Young. Effects of disease on clinical lab. tests 3 rd ed. AACC Pres. 1997

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Bilimsel Tıbbi Ürünler Paz. San. Tic. Ltd. Şti.