REF	CONT	REF	CONT	REF	CONT
IRO-10150	3x40mL+1x30mL	IRO-10228M	4x45 mL+2x24mL		
IRO-10400	4x80mL+1x80mL	IRO-10228P	4x45 mL+2x24mL		
IRO-10160A	4x32mL+4x8mL	IRO-10200M2	4x40 mL+2x20mL		
IRO-10102A2	4x20mL+2x11mL	IRO-10200M3	4x40 mL+2x20mL		

: 2-8 °C

**KULLANIM AMACI**

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
INTENDED USE

For the quantitative determination of iron in serum.

METHOD HISTORY

Iron measurements are used in the diagnosis and treatment of diseases and conditions such as iron deficiency anemia, hemochromatosis (characterized by a progressive increase in iron stores leading to organ impairment), chronic inflammatory disorders, hepatitis, and lead poisoning. Iron exists in the body as a component of hemoglobin and myoglobin as well as bound to transferrin for the transport in plasma and stored in ferritin.

Principle:

In an acidic medium transferrin bound iron dissociates into ferric ions. Reduced iron reacts with the Ferene to form a blue complex which absorbs at 595 nm. The absorbance is directly proportional to the serum iron concentration.

REAGENT COMPOSITION

Active Ingredients Concentration

Reagent 1	Concentration
Acetat Buffer	120 mmol/L
Thiourea	90 mmol/L
Surfactants	
pH 4.5 ± 0.1	
Reagent 2	Concentration
Ferene	3 mmol/L
Thiourea	20 mmol/L
Ascorbic Asid	

Precautions

- All reagents are toxic. Do not pipette by mouth. Avoid all contact.
- This reagent is for *in vitro* diagnostic use only.
- Use only disposable material to avoid iron contamination. Rinse glass material with diluted HCl and copious dist. water.

Reagent 1:

H317 May cause an allergic skin reaction

P302+352: IF ON SKIN: Wash with soap and water.

REAGENT STORAGE

Store all reagents refrigerated at 2-8°C.

REAGENT DETERIORATION

All reagents should be clear. Turbidity may indicate contamination and the reagent should not be used.

SPECIMEN COLLECTION AND STORAGE

- Fresh, unhemolyzed serum is the specimen of choice.
- Serum should be separated as soon as clot has formed.
- Heparinized plasma may be used but other anticoagulants should not be used to avoid possible iron contamination.
- Serum iron is reported to be stable for four days at room temperature (15-30°C) and seven days at 2-8°C.

INTERFERENCES

Bilirubin: No interference up to 10 mg/dL.

Lipemia: No interference in the presence of triglycerides up to 1000 mg/dL.

Ascorbic acid: No interference up to 20 mg/dL.

MATERIALS REQUIRED BUT NOT PROVIDED

- Accurate pipetting devices
- Test tubes/rack
- Timer
- Spectrophotometer able to read at 595 nm.
- Iron-free deionized water.
- Heating bath/block (37°C).

PROCEDURE (MANUAL)

Wavelength	: 595 nm
Working temperature	: 37°C
Optical path	: 1 cm

Assay type : Endpoint
Direction : Increasing

SERUM IRON

REAGENT 1	Blank	Calibrator	Sample
Distilled water	1000 µL	1000 µL	1000 µL
Iron Standard	100 µL	--	--
Sample	--	100 µL	--
	--	--	100 µL

- Mix, read the absorbance against blank in 1 - 5 minutes (Abs.1)

IRON CROMOGEN	Blank	Calibrator	Sample
	250 µL	250 µL	250 µL

- Mix and place all tubes in heating bath at 37°C for 10 minutes.

- After 10 minutes, read absorbance against blank.(Abs.2)

- ΔAbs.= (Abs.2 – Abs.1)

CALCULATION

Abs= Absorbance

$\frac{\Delta \text{Abs. Sample}}{\Delta \text{Abs. Calibrator}} \times \text{Conc. of Calibrator} = \text{Total Iron} (\mu\text{g/dL})$

ΔAbs. Calibrator

CALIBRATION

The assay requires the use of an Iron Standard or **bt products** Chem-Calibrator

QUALITY CONTROL

Serum controls with known normal and abnormal values should be run routinely to monitor the validity of the reaction.

EXPECTED VALUES

Total Iron

Women 37 – 165 µg/dL

Man 40 – 155 µg/dL

Children 35 – 135 µg/dL

It is strongly recommended that each laboratory determine the normal range for its particular population.

PERFORMANCE

Sensitivity/ Limit of Detection(LOD):

The lower limit of detection is 5,0 µg/d.

Linearity:

When the recommendation is designed and studied, the test is linear up to 1000 µg/dL. Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

Precision :

	Within run (n=20)		Between run (n=20)	
Mean (µg/dL)	102,4	230	112,8	257,5
SD	3,42	3,84	4,41	4,97
CV(%)	3,34	1,67	3,91	1,93

Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x).The results obtained using 40 samples were the following:

Correlation Coefficient (r) =0,9960

Regression $y = 0,9649 x + 3,61$

The results of the performance characteristics depend on the analyzer used.

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BİLİMSEL TIBBİ ÜRÜNLER PAZ.SAN.VE TİC.LTD.ŞTİ