
REF	CONT	REF	CONT	REF	CONT
LDH-30100	2x40mL+2x10mL	LDH-30228M	4x45mL+2x24mL	LDH-30100A2	4x20 mL+2x10mL
LDH-30250	5x40mL+1x50 mL	LDH-30228P	4x45mL+2x24mL		
LDH-30500	5x80mL+1x100 mL	LDH-30200M2	4x40mL+2x20mL		
LDH-30160A	4x32mL+4x8mL	LDH-30200M3	4x40mL+2x20mL		

KULLANIM AMACI

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: 2-8 °C

INTENDED USE

For the *in vitro* quantitative determination of lactate dehydrogenase activity in human serum.


CLINICAL SIGNIFICANCE

LDH is an enzyme that is found in highest concentration in the liver, kidney, heart, erythrocytes and skeletal muscle. LDH measurements are used in the diagnosis and treatment of liver diseases such as acute viral hepatitis, cirrhosis, and metastatic carcinoma of the liver, cardiac diseases such as myocardial infarction and tumors of the lung or kidneys.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

Principle:

Lactate dehydrogenase (LDH) catalyses the reduction of pyruvate by NADH, according to the following reaction:

The rate of decrease in concentration of NADPH, measured photometrically, is proportional to the catalytic concentration of LDH present in the sample¹.

REAGENT COMPOSITION

Active Ingredients

Concentrations

REAGENT 1

TRIS Buffer 100.0 mM
Pyruvate 3.0 mM
pH 7.8±0.2

REAGENT 2

NADH 1.6 mM
pH 10.0±0.3

Concentrations are those in the working reagent.

Precautions

- For *in vitro* diagnostic use only.
- DO NOT pipette by mouth. Avoid contact with skin and eyes. If spilled, thoroughly wash affected area with water.
- Reagent contains Sodium Azide as a preservative.
- Do not use the reagent after the expiration date printed on the kit.

Reagent 1:

H360 d May damage fertility or the unborn child
P308+313 If exposed or concerned: get medical advice / attention

REAGENT PREPARATION

Reagents are supplied in a two vial, ready to use, liquid form.
Working solution by mixing 4 parts of Reagent 1 with 1 part of Reagent 2 (e.g., 20mL R 1 to 5mL R 2).

REAGENT STORAGE

- Store the reagents refrigerated at 2–8°C.
- The reagents are stable until the expiration date, when stored at 2–8°C.
- Working reagent is stable for 2 weeks, when stored at 2–8°C.
- Do not freeze the reagents.

REAGENT DETERIORATION

The reagent should not be used if:

- The initial absorbance of the reagent is less than 1.00 when read against water at 340 nm.
- The reagent fails to meet stated parameters of performance.

SPECIMEN COLLECTION AND STORAGE

- Serum or heparinized plazma. Do not use oxalates as anticoagulants since they inhibit the enzyme.
- Non-hemolyzed serum is recommended. Red cells contain large concentrations of LDH.
- The serum should be removed from the clot promptly.
- Samples should be assayed soon after collection. LDH in serum is reported stable for 2 to 3 days at room temperature
- Do not freeze or expose the serum to high temperatures (37°C).

INTERFERENCES

Bilirubin: No interference up to 40 mg/dL.

Lipemia: No interference in the presence of triglycerides up to 1500 mg/dL.

Ascorbic acid: No interference up to 25 mg/dL.

Certain drugs and substances affect LDH activity.

ADDITIONAL EQUIPMENT REQUIRED BUT NOT PROVIDED

- A clinical chemistry analyzer capable maintaining constant temperature (37°C), and measuring absorbance at 340nm.
- Deionized water and related equipment, e.g.: pipettes
- Analyzer specific consumables, e.g.: sample and read cups
- Control materials

PROCEDURE

Wavelength : 340 nm
Temperature : 37°C
Optical path : 1 cm
Assay type : Kinetic
Direction : Decreasing

- Bring reagent to room temperature (15-30°C).
- Set the photometer to 0 (zero) absorbance with distilled water.

	Sample
Reagent 1	800 µL
Reagent 2	200 µL
Sample	20 µL

- Mix and incubate at 37° C.
- After 60 sec., read and record absorbance.
- Return tube to 37° C.
- Repeat readings every minute for the next three minutes
- Calculate the average absorbance difference per minute. (ΔAbs/Min.)

The ΔAbs/Min multiplied by the factor 8095 will yield results in IU/L.

Procedure Notes:

- The reagent and sample volumes may be altered proportionally to accommodate various instrument requirements.
- If the spectrophotometer being used is equipped with a temperature controlled cuvette, the reaction mixture may be left in the cuvette while the absorbance readings are taken.

CALCULATION

One Unit (U/L) is defined as the amount of enzyme that catalyzes the transformation of one micromole of substrate per minute.

U/L = ΔAbs/Min.x 8095

NOTE:

- If test parameters are altered, the factor has to be recalculated using the above formula.
- To convert to SI Units (nkat/L) multiply U/L by 16,67

CALIBRATION

The procedure is standardized by means of the millimolar absorptivity of NADH, taken as 6,22 at 340 nm under the test conditions described.

QUALITY CONTROL

The integrity of the reaction should be monitored by use of a two level control with LDH-P values.

EXPECTED VALUES

25 °C 120 – 240 U/L
30 °C 160 – 320 U/L
37 °C 230 – 460 U/L

It is strongly recommended that each laboratory determine the normal range for its particular population.

PERFORMANCE

Sensitivity/ Limit of Detection (LOD):

The lower limit of detection is 5 U/L.

Linearity:

When the recommendation is designed and studied, the test is linear up to 1500 U/L. Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

Precision :

	Within run (n=20)		Between run (n=20)	
Mean (U/L)	284,5	572,8	286,1	563,9
SD	5,76	6,14	5,40	6,54
CV(%)	2,02	1,07	1,89	1,16

Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x).The results obtained using 40 samples were the following:

Correlation Coefficient (r) =0,9983

Regression y =0,9839 x +4,0828

The results of the performance characteristics depend on the analyzer used.

REFERENCES

- Pes ce A. Lactate dehydrogenase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1124-117. 438.
- Young DS. Etfeds of drugs on Clinical Lab. Tests. 4th ed AACC Press. 1995.
- Young DS. Etfeds of disease on Clinical Lab. Tests. 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry. 3rd ed AACC 1999.
- Tietz N W et ai. Clinical Guide to Laboratory Tests. 3rd ed AACC 1995.

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BİLİMSEL TIBBİ ÜRÜNLER PAZ.SAN.VE TİC.LTD.ŞTİ