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(mg/dL)

INTENDED USE

This reagent is intended for the *in vitro* quantitative determination of low-density lipoprotein cholesterol (LDL-C) in human serum.

METHODOLOGY

When a measurement system where cholesterol esterase and cholesterol oxidase are used, this method separates LDL from other lipoprotein (HDL, VLDL and chylomicron) by inhibiting their reaction through the selective inhibitory function of inorganic and organic phosphorous compound as well as surface acting agent. This enables an exclusive and direct determination of the LDL cholesterol level.

REAGENT COMPOSITION

Reagent 1:

Detergent
Cholesterol esterase : ≤ 200.000 U/L
Cholesterol oxidase : ≤ 200.000 U/L
Peroxidase : ≤ 200.000 U/L
4-aminoantipyrine

Reagent 2:

Detergent
TOOS
Tris Buffer

REAGENT PREPARATION

Reagent 1 and Reagent 2 are ready to use.

Precautions and Warnings:

- For *in vitro* diagnostic use only.
- Do not pipette by mouth.
- Avoid contact with skin and eyes. If spilt, thoroughly, wash affected areas with water.
- All specimens used in the test should be considered potentially infectious.
- Do not use the reagent after the expiration date printed on the kit.

REAGENT STORAGE

- Store the reagent at 2-8°C.
- The reagent is stable until the expiration date printed on the label when stored at 2-8°C protected from light and contamination prevented during their use.

REAGENT DETERIORATION:

Do Not Use Reagent If:

- The reagent develops turbidity.
- The reagent fails to recover control values within the assigned range.

SPECIMEN COLLECTION AND HANDLING

Use fresh patient serum or plasma samples (EDTA, Citrate). Fasting and non-fasting samples can be used.

Storage: If not analyzed promptly, specimens may be stored at 2-8°C for up to 7 days.

INTERFERENCES

Bilirubin: No interference up to 20 mg/dL.
Hemoglobin: No interference up to 500 mg/dL.
Lipemia: No interference in the presence of triglycerides up to 1000 mg/dL.
Ascorbic acid: No interference up to 25 mg/dL.

MATERIALS REQUIRED BUT NOT PROVIDED

- A clinical chemistry analyzer or photometer.
- Analyzer specific consumables, eg: sample cups.
- Control and calibrator materials.

PROCEDURE

Primary Wavelength : 600 nm
Secondary Wavelength : 700 nm
Working temperature : 37°C
Optical path : 1 cm
Assay type : Endpoint/sample blank
Direction : Increasing

- Bring reagent to room temperature (15-30°C).

	Standard	Sample
Reagent 1	300 µL	300 µL
Standard	4 µL	--
Sample	--	4 µL

- Mix and then incubate for 5 min at 37°C, read the absorbance of the sample (Abs.1 sample) and standard (Abs.1 standard) against distilled water.

	Standard	Sample
Reagent 2	100 µL	100 µL

- Return tube to 37°C and after 5 min. read the final absorbance of the sample (Abs.2 sample) and standard (Abs.2 standard) against distilled water.

- Determine the Δ Abs. for Sample and Standard.

CALCULATIONS

(ΔAbs. = Abs.2 - Abs.1)

$\Delta \text{Abs. Sample} \times \text{Concentration of standard} = \text{LDL-Cholesterol (mg/dL)}$

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Limitations

- Anticoagulants containing heparin should not be used.
- Protect the reagents from direct sunlight.
- Store the reagents as per instructions.

CALIBRATION

bt products Chem-Calibrator or HDL/LDL Calibrator is required for calibration. If control results are found to be out of range, the procedure should be recalibrated.

QUALITY CONTROL

Reliability of test results should be routinely monitored with control materials that reasonably emulate performance of patient specimens. Quality control materials are intended for use only as monitors of accuracy and precision. Controls should be run with every working shift in which LDL-Cholesterol assays are performed.

RESULTS

LDL-Cholesterol concentration is expressed as mg/dL.

To convert from conventional units to S.I.units, multiply the conventional units by 0.02586.

$\text{mg/dL} \times 0.02586 = \text{mmol/L}$ LDL-Cholesterol

$\text{mmol/L} \times 38.66 = \text{mg/dL}$

EXPECTED VALUES

<100 mg/dL Optimal
100 - 129 mg/dL Near optimal / above optimal
130 - 159 mg/dL Borderline high
160 - 189 mg/dL High
≥190 mg/dL Verv high

Each laboratory must establish its own range of expected values.

PERFORMANCE

Sensitivity/ Limit of Detection(LOD):

The lower limit of detection is 0,20 mg/dL.

Linearity:

When the recommendation is designed and studied, the test is linear up to 250 mg/dL. Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

Precision :

	Within run (n=20)		Between run (n=20)	
Mean (mg/dL)	41,8	78,2	41,3	76,0
SD	0,81	1,98	0,78	1,50
CV(%)	1,95	2,54	1,89	1,97

Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x).The results obtained using 40 samples were the following:

Correlation Coefficient (r) =0,9953

Regression $y=0,9199x - 2,9733$

The results of the performance characteristics depend on the analyzer used.

REFERENCES

- "Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult treatment Panel III)", JAMA, 285:2486 (2001).
- Crouse JR et al., Studies of low density lipoprotein molecular weight in human beings with coronary artery disease, J. Lipid Res., 26: 566 (1985).
- Badimon 11, Badimon L., Fuester V., Regression of Atherosclerotic Lesions by High-Density Lipoprotein Plasma Fraction in the Cholesterol-Fed Rabbit, Journal of Clinical Investigation, 85: i 234-4 i (1990).
- Castelli WP et al., Cholesterol and other lipids in coronary heart disease, Circulation, 55: 767 (1977).
- Barr DP, Rillss EM, Eder HA, Protein-lipid relationships in human plasma, Am. J. Med., 1 1:480 (1951).
- Gordon T. et al., High density lipoprotein as a protective factor against coronary heart disease, Am.J.Med.,62:707 (1977).
- William P., Robinson D., Baily A., High density lipoprotein and coronary risk factor, Lancet, 1;72 (1979).
- Kannel WB, Castelli WP, Gordon T., Cholesterol in the prediction of atherosclerotic disease; New perspectives based on the Framingham study, Am. Intern. Med.90;85 (1979).
- Hongbing Xiao Method and composition for determining low density lipoprotein cholesterol, Chinese Patent CN1379234A (2002).
- Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of the High Blood Cholesterol in Adults (Adult Treatment Panel III).
- Pisani T., Gebski CP, Leary Et, et al. Accurate Direct Determination of Low-Density Lipoprotein Cholesterol Assay. Arch Pathol Lab Med 1995; 119:1127.