REF	CONT	REF	CONT
MAG-50100	2 x 50 mL	MAG-50180M	4 x 45 mL
MAG-50300	6 x 50 mL	MAG-50180P	4 x 45 mL
MAG-50600	6 x 100 mL	MAG-50240M2	6 x 40 mL
MAG-50040A	4 x 10 mL	MAG-50240M3	6 x 40 mL

REF	CONT
MAG-50100A2	5 x 20 mL

**KULLANIM AMACI**

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✂ : 2-25 °C


INTENDED USE

For the quantitative *in vitro* determination of magnesium in serum. A colorimetric xylidyl blue complex method

Principle

Magnesium ion forms a colored complex with xylidyl blue under alkaline conditions. The intensity of the developed color is proportional to the magnesium ion concentration of the sample.

REAGENT COMPOSITION

Active ingredients

Active ingredients	Concentration
Tris Buffer	200 mM
Xylidyl blue	0,1 mM
Detergent	

Precautions:

- This reagent is for *in vitro* diagnostic use only.
- Do not pipette by mouth. Avoid contact with skin and clothing.
H314 Causes severe skin burns and serious eye damage.
P303+P361+P353 IF ON SKIN: Remove/take off immediately all contaminated clothing. Rinse skin with water/shower
P304+P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
P305+P351+P338 IF IN EYES: Rinse continuously with water for several minutes.

Remove contact lenses if present and easy to do. Continue rinsing. Wear suitable protective equipment.

REAGENT PREPARATION

The reagent is ready for use.

REAGENT STORAGE

- Store reagent between 2-25°C.
- The reagent is stable until the expiration date stated on the label.

REAGENT DETERIORATION

Do not use reagent if:

- Reagent fails to achieve established values of fresh control sera.
- Reagent becomes visibly turbid.
- Reagent has an absorbance greater than 1,2 when measured against water at 500(480-520) nm.

SPECIMEN COLLECTION AND STORAGE

- Use fresh, unhemolyzed serum.
- Red cells contain twice the magnesium concentration as serum. A hemolyzed sample would falsely elevate results.
- Grossly icteric or lipemic specimens should not be used in this method.
- Urine (diluted in ratio of 1:10 with distilled water).

The sample should be adjusted to pH 3 - 4 with diluted hydrochloric acid.

INTERFERENCES

Bilirubin: No interference up to 10 mg/dL.
Hemoglobin: No interference up to 500 mg/dL.
Lipemia: No interference in the presence of triglycerides up to 1500 mg/dL.
Ascorbic acid: No interference up to 25 mg/dL.
A number of drugs and substances affect the concentration of magnesium.

ADDITIONAL EQUIPMENT REQUIRED BUT NOT PROVIDED

- A clinical chemistry analyzer capable maintaining constant temperature (37°C), and measuring absorbance at 500 nm (480-520).
- Deionized water and related equipment, e.g.: pipettes
- Analyzer specific consumables, e.g.: sample cups
- Control and calibrator materials.

PROCEDURE

System Parameters

Magnesium	
Temperature	: 37°C
Wavelength	: 500 nm(480 – 520)
Optical Path	: 1 cm
Assay Type	: Endpoint
Direction	: Increase

	Blank	Standard	Sample
Reagent	1000 µL	1000 µL	1000 µL
Distilled water	10 µL	--	--
Standard	--	10 µL	--
Sample	--	--	10µL

Mix and then incubate for 5 min. at 37°C. Measure the absorbance of sample and standard against the blank.

Procedure Note:

- The dilution factor should be considered at the calculation of magnesium ion concentration of urine.
- Do not use reagents after the expiry date stated on each reagent container label.
- Do not use products, test solutions and reagents described above for any purpose other than described herein.

CALCULATIONS:

(Abs= Absorbance)

$\frac{\text{Abs. Sample}}{\text{Abs. Standard}} \times \text{Concentration of std} = \text{Magnesium (mg/dL)}$

CALIBRATION

Use an aqueous **bt products** Magnesium standard or an appropriate serum calibrator.

QUALITY CONTROL

The integrity of the reaction should be monitored by use of a two level control with known Magnesium values

EXPECTED VALUES

Serum: 1.6 – 2.5 mg/dL (0.65-1.03 mmol/L)

It is highly recommended that each laboratory establish its own reference range.

NOTE: "mg/dL" may be converted to "mEq/L" by dividing the result by 1,21525.

PERFORMANCE

Sensitivity/ Limit of Detection(LOD):

The lower limit of detection is 0,16 mg/dL.

Linearity:

When the recommendation is designed and studied, the test is linear up to 5 mg/dL. Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

Precision :

	Within run (n=20)	Between run (n=20)		
Mean (mg/dL)	1,93	2,65	1,92	2,55
SD	0,03	0,04	0,03	0,04
CV(%)	1,76	1,50	1,35	1,76

Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x).The results obtained using 40 samples were the following:

Correlation Coefficient (r) =0,9969

Regression $y=1,1124x+0,0363$

The results of the performance characteristics depend on the analyzer used.

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BİLİMSEL TIBBİ ÜRÜNLER PAZ.SAN.VE TİC.LTD.ŞTİ