REF	CONT	REF	CONT	REF	CONT
PHO-10100	2 x 50 mL	PHO-10180M	4 x 45 mL	PHO-10100A2	5 x 20 mL
PHO-10300	6 x 50 mL	PHO-10180P	4 x 45 mL		
PHO-10600	6 x 100 mL	PHO-10240M2	6 x 40 mL		
PHO-10040A	4 x 10 mL	PHO-10240M3	6 x 40 mL		

KULLANIM AMACI

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INTENDED USE

For the *in vitro* quantitative determination of Inorganic Phosphorus in serum.

METHODOLOGY

The measurement of Inorganic phosphorus in serum is usually accomplished by forming a phosphomolybdate complex and in turn reducing it to a molybdenum blue color complex. Methods differ as to the choice of reducing agents: stannous chloride, phenylhydrazine, aminonaphtholsulfonic acid, ascorbic acid, p-methylamino phenolsulfate, N-phenyl-p-phenylenediamine and ferrous sulfate. These methods suffered from color instability, deproteinization steps and complexity of performance. The addition of a surfactant eliminated the need to prepare a protein-free filtrate, accelerated color production, stabilized the color and simplified the procedure. Many of the components in these reagents were unstable and had to be stored separately. The quantitative measurement of unreduced phosphomolybdate complexes was first reported by Simonsen in 1946. Daly and Ertingshausen adapted that technique for the determination of Inorganic phosphorus in 1972. Amador and Urban modified this procedure further the same year. The present method is a modification of the above procedure using a single, stable reagent performing in the UV range.

Principle

Inorganic Phosphorus + H₂SO₄ + Ammonium Molybdate → Unreduced Phosphomolybdate Complex.

Inorganic phosphorus reacts with ammonium molybdate in an acid medium to form a phosphomolybdate complex which absorbs light at 340 nm. The absorbance at this wavelength is directly proportional to the amount of inorganic phosphorus present in the sample.

REAGENT COMPOSITION

Active Ingredients	Concentration
Ammonium Molybdate	0.81 mM
Sulfuric Acid	255 mM
Surfactant.	
pH 0.4 ± 0.2	

Precautions

- This reagent is for *in vitro* diagnostic use only.
- This reagent is an acid and is caustic, Avoid contact with skin. Flush with plenty of water if contact occurs, do not pipette by mouth.

Reagent contains Sulfuric acid and Triton x-100.

H314 Causes severe skin burns and serious eye damage.

P303+P361+P353 IF ON SKIN: Remove/take off immediately all contaminated clothing. Rinse skin with water/shower

P304+P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.

P305+P351+P338 IF IN EYES: Rinse continuously with water for several minutes

.Remove contact lenses if present and easy to do. Continue rinsing

Wear suitable protective clothing, gloves and eye/face protection.

REAGENT PREPARATION

Reagent is ready to use.

REAGENT STORAGE

Store the reagent at 2-25°C.

REAGENT DETERIORATION

Do not use the reagent if:

- Reagent read against water has an absorbance greater than 0,500 at 340 nm.
- The reagent fails to meet stated parameters of performance.

SPECIMEN COLLECTION AND STORAGE

- Unhemolyzed serum is specimen of choice, serum should be removed from the red cell clot as soon as possible.
- Plasma should not be used since anticoagulants may produce falsely low values.
- Serum inorganic phosphorus is stable for one week refrigerated and for three weeks frozen.

INTERFERENCES

Bilirubin: No interference up to 40 mg/dL.

Lipemia: No interference in the presence of triglycerides up to 1500 mg/dL.

Ascorbic acid: No interference up to 25 mg/dL.

A number of drugs and substances may affect the accuracy of phosphorus.

ADDITIONAL EQUIPMENT REQUIRED BUT NOT PROVIDED

- A clinical chemistry analyzer capable maintaining constant temperature (37°C), and measuring absorbance at 340nm.
- Deionized water and related equipment, e.g.: pipettes
- Analyzer specific consumables, e.g.: sample cups
- Control and calibrator materials.

PROCEDURE

Wavelength : 340 nm
Working temperature : 37°C

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Optical path : 1 cm
Assay type : Endpoint
Direction : Increasing

Reagent	Blank	Standard	Sample
1000 µL	1000 µL	1000 µL	1000 µL
Distilled water	10 µL	--	--
Standard	--	10 µL	--
Sample	--	--	10 µL

Mix and then incubate for 5 min. at 37°C. Measure the absorbance of sample and standard against the reagent blank.

Procedure Note:

The reagent and sample volumes may be altered proportionally to accommodate various instrument requirements.

CALCULATIONS:

(Abs = Absorbance)

$\frac{\text{Abs. Sample}}{\text{Abs. Standard}} \times \text{Concentration of standard (mg/dL)} = \text{Phosphorus (mg/dL)}$

LIMITATIONS:

- Samples with values exceeding 15 mg/dL should be diluted 1:1 with saline and re-run. The final answer should be multiplied by two.
- Detergents containing phosphate should not be used for cleaning glassware used in this procedure.

CALIBRATION

Use an aqueous Phosphorus standard or an appropriate serum calibrator.

QUALITY CONTROL

The integrity of the reaction should be monitored by use of a two level control with known Phosphorus values

EXPECTED VALUES

Adults : 2.5 - 4.8 mg/dL

Children: 4.0 - 7.0 mg/dL

Values are decreased during menstrual period and after meals.

It is strongly recommended that each laboratory establish its own normal range.

PERFORMANCE

Sensitivity/ Limit of Detection(LOD):

The lower limit of detection is 0,077 mg/dL.

Linearity:

When the recommendation is designed and studied, the test is linear up to 15 mg/dL. Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

Precision :

	Within run (n=20)		Between run (n=20)	
Mean (mg/dL)	4,13	7,99	4,11	8,12
SD	0,07	0,14	0,08	0,21
CV(%)	1,67	1,70	1,93	2,65

Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x). The results obtained using 40 samples were the following:

Correlation Coefficient (r) =0,9960

Regression $y = 0,8355x + 0,3051$

The results of the performance characteristics depend on the analyzer used.

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BİLİMSEL TIBBİ ÜRÜNLER PAZ.SAN.VE TİC.LTD.ŞTİ