POT-10120  
POT-10060A  
POT-10050A23x30mL+1x30mL+2x2mL Std.  
2x22,5mL+2x7,5mL+2x2mL Std.  
3x12,5mL+1x12,5mL+2x2mL Std.POT-10080M2  
POT-10080M32x30mL+2x10mL+2x2mL Std.  
2x30mL+2x10mL+2x2mL Std.**KULLANIM AMACI**

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**BİLİMSEL KAYNAKÇA:**

- M.N. Berry,R. D. Mazzachi, M. Pejakovic, and M. J. Peake Enzymatic Determination of Potassium in Serum. *CLIN. CHEM* 35/5, 817-820
- Wu, A.H.B.ed.Tietz clinical guide to laboratory tests, 41h edition, p. 880. W.B. Saunders Company, St. Louis (2006).
- Bergmeyer,-H.U.,Gawehn, K., and Grassl, M. (1974) in *Methods of Enzymatic Analysis*. Second Edition, Volume I, 509-510, Academic Press, Inc., New York.

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⚠ : 2-8 °C

### INTENDED USE

For the quantitative determination of Potassium in human serum and plasma.

### METHODOLOGY

Potassium is determined spectrophotometrically through a kinetic coupling assay system using potassium dependent pyruvate kinase. Pyruvate generated is converted to lactate accompanying conversion of NADH analog to NAD analog. The corresponding decrease of optical density at 380nm is proportional to the potassium concentration in the serum.

### REAGENT COMPOSITION

#### Reagent 1

LDH ≥ 50 U/mL  
NADH ≥ 0,15 mM

Buffer, Azide and stabilizers.

#### Reagent 2

Pyruvate Kinase ≥ 0,9 mM

Azide and stabilizers.

### Standards

SPOT1-102 1x2 mL  
SPOT2-102 1x2 mL

### Precautions:

1. This reagent is for *in vitro* diagnostic use only.
2. Do Not Ingest. Avoid contact with skin and eyes.
3. These reagents contain lithium azide. Azide compounds may react with lead or copper plumbing to form potentially explosive compounds. Flush drains with CO

### REAGENT PREPARATION

All the reagents are ready to use.

### REAGENT STORAGE

Reagent 1 and Reagent 2 are provided in ready to use liquid form and are stable up to expiration date when stored under 2-8°C.

Do not use reagents over the expiration date.

### REAGENT DETERIORATION

Do not use the reagent if:

1. Presence of particles and turbidity.

### SPECIMEN COLLECTION AND STORAGE

For use non-hemolysed serum .

Stability: Serum samples should be collected such that testing is performed as soon as possible and within 5 days after the specimen collection.

### INTERFERANCES

The assay is not interfered by the following substances at indicated concentrations: Na<sup>+</sup> 150 mM, NH<sub>4</sub><sup>+</sup> 0,5 mM, Ca<sup>2+</sup> 7,5 mM, Pi 2.0 mM, ascorbic acid 10,0 mM, Zn<sup>2+</sup> 0,5 mM, Fe<sup>3+</sup> 0,5 mM, Cu<sup>2+</sup> 0,5 mM, triglycerides 1000 mg/dL, hemoglobin 500 mg/dL, conjugated bilirubin 20 mg/dL and unconjugated bilirubin 15 mg/dL.

### Materials Required but not Provided

1. Spectrophotometer or colorimeter measuring at 380-405 nm.
2. Matched cuvettes 1.0 cm light path.
3. General laboratory equipment.

### PROCEDURE

Wavelength : 380 nm  
Assay Type : Fixed time  
Reaction Direction : Decreasing  
Temperature : 37°C  
Optical path : 1 cm

1. Bring to room temperature ( 15 -30 °C)
2. Set the photometer to 0 (zero ) with distilled water.

	Sample or Standard tube
R1	225 µL
Sample or Standard	6 µL

3. Mix and incubate for 5 minutes at 37° C

R2	75 µL
----	-------

4. Mix and incubate for 1 minutes at 37° C and read the absorbance immediately (Abs.1)

5. Return tube to 37 °C and after 3 minutes read the final Abs.(Abs.2)

6. Calculation  $\Delta A = A1 - A2$ :

### CALIBRATION

This assay should be calibrated using low and high potassium standards. Potassium concentration in sample is determined from linear calibration curve using SPOT1-102 and SPOT2-102 potassium standards.

For analyzers that require a zero standart, 0,9 g/dL NaCl (saline) can be used for that purpose and SPOT1-102 and SPOT2-102 potassium standards are used as standard 2 and standard 3 respectively. Weekly calibration is recommended

### CALCULATIONS

Calculate the absorbance difference (Abs.1 – Abs.2) of each point of the calibration curve and plot the values obtained against the Potassium concentration of each standard dilution. Potassium concentration in the sample is calculated by interpolation of its (Abs.1 – Abs.2) in the calibration curve

### QUALITY CONTROL

The validity of the reaction should be monitored using control sera with known normal and abnormal values. If control values are found outside the defined range check the instrument, reagent and calibration for problems.

Each laboratory should establish its own Quality control scheme and corrective actions if controls do not meet the acceptable tolerances.

### REFERENCE VALUES

3,5 - 5,1 mmol/L ( 13,7 – 19,9 mg/dL)

It is recommended that each laboratory establish its own reference range.

### PERFORMANCE

1. **Linearity:** When run as recommended the assay is linear from the measuring range of 2,0 mmol/L- 8.0 mmol/L
2. **Sensitivity:** The minimum detectable concentration of potassium with an acceptable level of precision was determined as 0,87 mmol/L.
3. **Precision:**

Within run:

	4,46 mmol/L K <sup>+</sup> 20 days, n=80)	6,86 mmol/L K <sup>+</sup> 20 days, n=80)
Mean(mmol/L)	4,62	6,96
CV(%)	1,12	1,20

Total :

	4,46 mmol/L K <sup>+</sup> 20 days, n=80)	6,86 mmol/L K <sup>+</sup> 20 days, n=80)
Mean(mmol/l)	4,62	6,96
CV(%)	1,77	1,77

### 4. Accuracy

The potassium assay was tested on the Autoanalyser and results obtained were compared to an ISE method. A total of 52 serum samples ranging from 2.7-7.7 mmol/L potassium and two sets of serum based controls were tested in both assays. The above described accuracy study showed that the BTPRODUCTS reagent method had good correlation with existing ISE method with a correlation coefficient of 0.98 with a slope of 1.07 and -0.30 intercept.

### REFERENCES

1. M.N. Berry, R. D. Mazzachi, M. Pejakovic, and M. J. Peake Enzymatic Determination of Potassium in Serum. *CLIN. CHEM* 35/5, 817-820
2. Wu, A.H.B.ed. Tietz clinical guide to laboratory tests, 41h edition, p. 880. W.B. Saunders Company, St. Louis (2006).
3. Bergmeyer, H.U., Gawehn, K., and Grassl, M. (1974) in *Methods of Enzymatic Analysis*. Second Edition, Volume I, 509-510, Academic Press, Inc., New York.

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Bilimsel Tıbbi Ürünler Paz.San.Tic.Ltd.Şti.