
REF	CONT	REF	CONT	REF	CONT
PRO-10100	2 x 50 mL	PRO-10180M	4 x 45 mL	PRO-10100A2	5 x 20 mL
PRO-10300	6 x 50 mL	PRO-10180P	4 x 45 mL		
PRO-10600	6 x 100 mL	PRO-10240M2	6 x 40 mL		
PRO-10080A	4 x 20 mL	PRO-10240M3	6 x 40 mL		

### KULLANIM AMACI

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### INTENDED USE

For the *in vitro* quantitative determination of total protein concentration in serum.

### METHODOLOGY

The color reaction of protein molecules with cupric ions, known as the Biuret color reaction, has been known since 1878. Since the Riegler publications of 1914, several attempts have been made to stabilize the cupric ions in the alkali reagent. Kingsley, modified the procedure in 1939 and 1942 to include the use of sodium potassium tartrate as a complexing agent. This procedure was later modified by Weichselbaum and Gornall. The present method is based on these modifications

**Principle:** Alkali Medium

Protein + Cu<sup>++</sup> → Colored Complex

Protein in serum forms a blue colored complex when reacted with cupric ions in an alkaline solution. The intensity of the violet color is proportional to the amount of protein present when compared to a solution with known protein concentration.

### REAGENT COMPOSITION

Active Ingredients	Concentrations
Sodium Hydroxide	600 mM
Cupric Sulfate	12 mM
Potassium Sodium Tartrate	32 mM
Potassium Iodide	30 mM

### Precautions

1. This reagent is for *in vitro* diagnostic use only.
2. Avoid ingestion. Do Not Pipette By Mouth. In case of ingestion drink large amounts of water and seek medical attention quickly.
3. Avoid contact with skin and eyes. The reagent contains sodium hydroxide which is corrosive. In case of contact with skin, flush with water. For eyes, seek medical attention.

Reagent contains Sodium hydroxide.

H314 Causes severe skin burns and serious eye damage.

P303+P361+P353 IF ON SKIN: Remove/take off immediately all contaminated clothing. Rinse skin with water/shower

P304+P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.

P305+P351+P338 IF IN EYES: Rinse continuously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing

Wear suitable protective clothing, gloves and eye/face protection.

### REAGENT PREPARATION AND DETERIORATION

Reagent is supplied ready to use. The reagent should be a clear pale blue solution. Turbidity or the presence of a black precipitate indicates reagent deterioration and should not be used.

### REAGENT STORAGE

Store reagent at room temperature. (2-25°C)

### SPECIMEN COLLECTION AND STORAGE

1. Unhemolyzed serum is the specimen of choice.
2. Gross hemolysis will cause elevated results because of the released hemoglobin as well as the increase in background color.
3. Lipemic sera cause elevated results and should be run with a serum blank.
4. Samples with bromosulfophthalein (BSP) will result in falsely elevated results.
5. Protein in serum is stable for one week at room temperature (18 - 25°C) and for at least one month refrigerated (2 - 8°C) when guarded against evaporation.

### INTERFERENCES

Bilirubin: No interference up to 40 mg/dL.

Hemoglobin: No interference up to 200 mg/dL.

Lipemia: No interference in the presence of triglycerides up to 1500 mg/dL.

Ascorbic acid: No interference up to 25 mg/dL.

### ADDITIONAL EQUIPMENT REQUIRED BUT NOT PROVIDED:

1. A clinical chemistry analyzer capable maintaining constant temperature (37°C), and measuring absorbance at 540nm.
2. Iron-free deionized water and related equipment, e.g.: pipettes
3. Analyzer specific consumables, e.g.: sample and read cups
4. Control and calibrator materials

### PROCEDURE

Wavelength : 540 nm

Working temperature : 37°C  
Optical path : 1 cm  
Assay type : Endpoint  
Direction : Increasing

	Blank	Standard	Sample
Reagent	1000 µL	1000 µL	1000 µL
Distilled water	20 µL	--	--
Standard	--	20 µL	--
Sample	--	--	20 µL

Mix and then incubate for 5 minutes at 37°C. Measure the absorbance of sample and standard against the reagent blank.

### Procedure Notes

1. Final color is stable for 60 minutes at room temperature.
2. The reagent and sample volumes may be altered proportionally to accommodate various instrument requirements.

### CALCULATION

(Abs = Absorbance)

Abs. of Sample x conc. of std. (g/dL) = Total protein (g/dL)

Abs. of Standard

### LIMITATIONS

The Biuret procedure is not sensitive at low ranges (<1g/dL). Do not use for urine or spinal fluid.

### CALIBRATION

Use an aqueous Total protein standard or an appropriate serum calibrator.

### QUALITY CONTROL

The integrity of the reaction should be monitored by use of a two level control with known Total Protein values.

### EXPECTED VALUES

6,2 – 8,5 g/dL

The effect of posture, when blood is drawn, varies with the individual but recumbent values are usually lower than ambulatory. Differences may be as much as 1,2 g/dL. It is strongly recommended that each laboratory establish its own normal range.

### PERFORMANCE

#### Sensitivity/ Limit of Detection(LOD):

The lower limit of detection is 0,061 g/dL.

#### Linearity:

When the recommendation is designed and studied, the test is linear up to 10 g/dL. Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

#### Precision :

	Within run (n=20)		Between run (n=20)	
Mean (g/dL)	4,55	7,55	4,51	7,68
SD	0,04	0,08	0,07	0,11
CV(%)	0,97	1,08	1,61	1,38

#### Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x). The results obtained using 40 samples were the following:

Correlation Coefficient (r) =0,9975

Regression equation y=1,0515 x +0,0225

The results of the performance characteristics depend on the analyzer used.

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BİLİMSEL TIBBİ ÜRÜNLER PAZ.SAN.VE TİC.LTD.ŞTİ