
REF	CONT	REF	CONT	REF	CONT
TRI-10100	2 x 50 mL	TRI-10300M	6 x 50 mL	TRI-10100A2	5 x 20 mL
TRI-10300	6 x 50 mL	TRI-10300P	6 x 50 mL		
TRI-10600	6 x 100 mL	TRI-10240M2	6 x 40 mL		
TRI-10180A	4 x 45 mL	TRI-10240M3	6 x 40 mL		

KULLANIM AMACI

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




INTENDED USE

For the *in vitro* quantitative determination of Triglycerides in serum.

METHODOLOGY

This method uses a modified Trinder color reaction to produce a fast, linear, endpoint reaction

Principle:

Triglycerides in the sample are hydrolyzed by lipase to glycerol and fatty acids. The glycerol is then phosphorylated by adenosine-5-triphosphate (ATP) to glycerol-3-phosphate (G₃P) and adenosine-5-diphosphate in a reaction catalyzed by glycerol kinase (GK). Glycerol-3-phosphate is then converted to dihydroxyacetone phosphate (DAP) and hydrogen peroxide by glycerophosphate oxidase (GPO). The hydrogen peroxide then reacts with 4-aminophenazone (4-AP) and p-chlorophenol in a reaction catalyzed by peroxidase to yield a red colored quinoneimine dye.

The intensity of the color produced is directly proportional to the concentration of Triglycerides in the sample when measured at 505 nm.

REAGENT COMPOSITION

Active Ingredients

Active Ingredients	Concentrations
ATP	0,1 mM
4-Aminophenazone	0,1 mM
p-chlorophenol	2,0 mM
GPO	3500 U/L
Lipoprotein Lipase	150000 U/L
GK	500 U/L
GOOD	50 mM

pH 6,3 ± 0,1

Precautions

1. Reagent contains Sodium Azide as a preservative.
2. This reagent is for *in vitro* diagnostic use only.

REAGENT PREPARATION

Reagent comes in a ready to use form.

REAGENT STORAGE

Store the reagent at 2-8°C (refrigerated).

The reagent is stable until the expiration date when stored at 2-8°C.

REAGENT DETERIORATION

Do not use the reagent if:

1. Reagent is turbid.
2. The reagent has an absorbance greater than 0,300 against water at 505 nm.
3. The reagent fails to recover stated values in control sera.

SPECIMEN COLLECTION AND STORAGE

1. Serum or plasma.
2. Triglycerides in serum is stable for 5 days at 2-8°C.

INTERFERENCES

Bilirubin: No interference up to 20 mg/dL.

Hemoglobin: No interference up to 200 mg/dL.

ADDITIONAL EQUIPMENT REQUIRED BUT NOT PROVIDED

1. A clinical chemistry analyzer capable maintaining constant temperature (37°C), and measuring absorbance at 490 – 550 nm.
2. Deionized water and related equipment, e.g.: pipettes
3. Analyzer specific consumables, e.g.: sample and read cups
4. Control and calibrator materials.

PROCEDURE

Wavelength	: 510 nm(490-550)
Working temperature	: 37°C
Optical path	: 1 cm
Assay type	: Endpoint
Direction	: Increasing

	Blank	Standard	Sample
Reagent	1000 µL	1000 µL	1000 µL
Distilled water	10 µL	--	--
Standard	--	10 µL	--
Sample	--	--	10µL

Mix and then incubate for 5 min at 37°C or 10 min.15-25 °C. Measure the absorbance of sample and standard against the reagent blank. The colour is stable at least 30 min.

PROCEDURE NOTES

1. The reagent and sample volumes may be altered proportionally to accommodate various instrument requirements.
2. Final color is stable for 30 minutes.

CALCULATION

Abs. = Absorbance

$\frac{\text{Abs. Sample}}{\text{Abs. Standard}} \times \text{Conc. Of standard} = \text{Triglycerides (mg/dL)}$

NOTE:To obtain values in S.I. units multiply mg/dL x 0.113 = mmol/L

LIMITATIONS

Glycerol (free glycerol and glycerol released upon hydrolysis of triglycerides) is measured by this procedure. Free glycerol levels in serum are generally low, but elevations may be caused by improper storage or sample contamination

CALIBRATION

Calibration with the standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum calibrator.

QUALITY CONTROL

The integrity of the reaction should be monitored by use of a two level control with known Triglycerides values.

EXPECTED VALUES

Men 40 – 160 mg/dL

Women 35 – 135 mg/dL

The expected range is taken from reference literature. Each laboratory should establish its own normal range.

PERFORMANCE

Sensitivity/ Limit of Detection(LOD):

The lower limit of detection is 0,63 mg/dL.

Linearity:

When the recommendation is designed and studied, the test is linear up to 1000 mg/dL. Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

Precision :

	Within run (n=20)		Between run (n=20)	
Mean (mg/dL)	125,6	215,3	121,6	208,9
SD	1,53	3,36	3,07	3,92
CV(%)	1,22	1,56	2,52	1,87

Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x).The results obtained using 40 samples were the following:

Correlation Coefficient (r) =0,9988

Regression equation y =0,9127x +1,69

The results of the performance characteristics depend on the analyzer used.

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BİLİMSEL TIBBİ ÜRÜNLER PAZ.SAN.VE TİC.LTD.ŞTİ