REF	CONT
UIBC-50150	2x50 mL+1x25 mL+1x25 mL
UIBC-50480	4x80 mL+1x80 mL+1x80 mL
UIBC-50120A	4x20 mL+2x10 mL+2x10 mL
UIBC-50120A2	4x20 mL+1x20 mL+1x20mL

REF	CONT
UIBC-50270M	3x45 mL+2x17 mL+2x17 mL
UIBC-50270P	3x45 mL+2x17 mL+2x17 mL
UIBC-50160M2	4x30 mL+1x20 mL+1x20 mL
UIBC-50240M3	4x40 mL+2x20 mL+2x20 mL

⚠ : 2-8 °C

KULLANIM AMACI

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⚠ : 2-8 °C

INTENDED USE

For the quantitative determination of iron and total iron-binding capacity in serum.

Principle:

Total Iron-Binding Capacity (TIBC): A known amount of ferrous ions are added to serum at an alkaline pH. The ferrous ions bind with transferrin at unsaturated iron-binding sites. The additional unbound ferrous ions are measured using the ferrozine reaction. The difference between the amount of ferrous ions added and the unbound ions measured is the unsaturated iron-binding capacity (UIBC). The TIBC is equal to the serum iron concentration plus the UIBC.

REAGENTS

UIBC/TIBC BUFFER REAGENT: Tris 500mM, pH 8.1 with surfactant, Sodium Azide 0.05% (w/v) as preservative.

UIBC/TIBC COLOR REAGENT(Cromogen) : Ferrozine in hydroxylamine hydrochloride.

UIBC/TIBC STANDARD (500 µg/dL): Ferrous chloride in hydroxylamine hydrochloride.

Precautions

- All reagents are toxic. Do not pipette by mouth. Avoid all contact.
- UIBC buffer contains sodium azide and may react with lead and copper plumbing to form highly explosive metal azides.
- On disposal, flush with a large volume of water to prevent azide accumulation.
- This reagent is for *in vitro* diagnostic use only.

CROMOGEN and STANDARD

H317 May cause an allergic skin reaction.

P302+352: IF ON SKIN: Wash with soap and water .

REAGENT STORAGE

Store all reagents refrigerated at 2-8°C.

REAGENT DETERIORATION

All reagents should be clear. Turbidity may indicate contamination and the reagent should not be used.

SPECIMEN COLLECTION AND STORAGE

- Fresh, unhemolyzed serum is the specimen of choice.
- Serum should be separated as soon as clot has formed.
- Heparinized plasma may be used but other anticoagulants should not be used to avoid possible iron contamination.
- Serum iron is reported to be stable for four days at room temperature (15-30°C) and seven days at 2-8°C.

INTERFERENCES

Lipemia: No interference in the presence of triglycerides up to 1500 mg/dL.

Ascorbic acid: No interference up to 25 mg/dL.

- Certain drugs and other substances are known to influence circulating iron levels.
- Iron contained in hemoglobin does not react in this method, therefore, slight hemolysis will not interfere.
- To make tubes, pipettes, etc. iron free, they must be washed with hot, dilute (1:2) hydrochloric or nitric acid, followed by several rinsings with iron-free deionized or distilled water.

MATERIALS REQUIRED BUT NOT PROVIDED

- Accurate pipetting devices
- Test tubes/rack
- Timer
- Spectrophotometer able to read at 560 nm.
- Iron-free deionized water.
- Heating bath/block (37°C).

PROCEDURE (AUTOMATED)

Refer to specific instrument application instructions.

PROCEDURE (MANUAL).

Wavelength	: 560 nm
Working temperature	: 37°C
Optical path	: 1 cm
Assay type	: Endpoint
Direction	: Increasing

UIBC (Unsaturated Iron-Binding Capacity)

	Blank	Standard	Sample
UIBC/TIBC Buffer Reagent	1000 µL	1000 µL	1000 µL
Distilled water	500 µL	250 µL	--
UIBC/TIBC Standard	--	250 µL	250 µL
Sample	--	--	250 µL

- Zero spectrophotometer at 560nm with the blank.
- Read and record absorbances of all tubes. (Abs.1 reading).

	Blank	Standard	Sample
UIBC/TIBC Color Reagent (Cromogen)	250 µL	250 µL	250 µL

3. Mix and place all tubes in heating bath at 37°C for 10 minutes.

4. Zero instrument at 560nm with reagent blank.

5. Read and record absorbances of all tubes. (Abs.2 reading).

UIBC CALCULATIONS

Abs.= Absorbance Std.= Standard

$$\text{Conc. Of Std.} = \frac{\text{Abs.2 Sample} - \text{Abs.1 Sample}}{\text{Abs.2 Std.} - \text{Abs.1 Std.}} \times \text{Conc. of std.} = \text{UIBC } \mu\text{L/dL}$$

Example:

Std. Conc. = 500 µL/dL

Abs.2 Sample – Abs.1 Sample = 0.10

Abs.2 Std.– Abs. 1 Std = 0.40

Therefore: $500 - \frac{0.10}{0.40} \times 500 = \text{UIBC } (\mu\text{L/dL})$

$500 - (0.25 \times 500) = 375 \mu\text{L/dL (UIBC)}$

NOTE: The difference between Abs.1 Sample and Abs.2 Sample Test may sometimes be very small due to a high degree of insaturation of transferrin with iron. The sample should be diluted 1:1 with iron-free water and re-assayed. The result is then multiplied by two.

CALCULATION

TIBC (Total Iron-Binding Capacity)

Iron Level + UIBC = TIBC (µL/dL)

SI Unit Conversion µL/dL x 0.179 = µmol/L

CALIBRATION

The procedure is calibrated with iron standard (500 µg/dL) included in each kit.

QUALITY CONTROL

Serum controls with known normal and abnormal values should be run routinely to monitor the validity of the reaction.

EXPECTED VALUES

TIBC = 250 - 400 µL/dL

Iron Saturation = 20-55%

It is strongly recommended that each laboratory determine the normal range for its particular population.

PERFORMANCE

Sensitivity/ Limit of Detection(LOD):

The lower limit of detection is 25 µL/dL.

Linearity:

When the recommendation is designed and studied, the test is linear up to 500 µL/dL. Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

Precision :

	Within run (n=20)		Between run (n=20)	
Mean (µL/dL)	168,1	209,1	172,7	216,4
SD	2,19	2,37	5,73	6,59
CV(%)	1,30	1,13	3,32	3,04

Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x).The results obtained using 40 samples were the following:

Correlation Coefficient (r) =0,9952

Regression $y = 1,0484 x - 4,761$

The results of the performance characteristics depend on the analyzer used.

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BİLİMSEL TIBBİ ÜRÜNLER PAZ.SAN.VE TİC.LTD.ŞTİ