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www.btproducts.com.tr - www.bilimseltip.com

REF	CONT	REF	CONT	REF	CONT
URE-11100	2x40mL+2x10mL	URE-11228M	4x45mL+2x24mL	URE-11100A2	4x20mL+2x10mL
URE-11250	5x40mL+1x50mL	URE-11228P	4x45mL+2x24mL		
URE-11500	5x80mL+1x100mL	URE-11300M2	6x40mL+3x20mL		
URE-11160A	4x32mL+4x8mL	URE-11200M3	4x40mL+2x20mL		

⚠ : 2-8 °C


INTENDED USE

For the quantitative determination of UREA in serum.

METHODOLOGY

Talke and Schubert introduced a totally enzymatic procedure in 1965 utilizing urease and glutamate dehydrogenase. The present procedure is based on a modification of their method.

Principle

Urea is hydrolyzed by urease to produce ammonia and carbon dioxide. The liberated ammonia reacts with α -ketoglutarate in the presence of NADH to yield glutamate. An equimolar quantity of NADH undergoes oxidation during the reaction resulting in a decrease in absorbance that is directly proportional to the urea concentration in the sample.

REAGENT COMPOSITION

Active Ingredients (R1+R2)

	Concentration
Buffer	
α -Ketoglutarate	< 4.0 mM
Urease	7500 U/L
GLDH	1000 U/L
ADP	> 2.0 mM
NADH	1,29 mM

Stabilizers and detergents

pH 7.4 \pm 0.1

REAGENT PREPARATION

Prepare working reagent by mixing 4 parts reagent 1 with 1 part reagent 2 (e.g. 200 μ L R1 with 50 μ L R2 reagent)

REAGENT STORAGE

1. Store reagents at 2-8°C.
2. Store working reagent at 2-8°C.
3. Working reagent is stable one week at 20-25°C and 4 weeks at 2-8°C.
4. Do not freeze the reagents

REAGENT DETERIORATION

The reagent should not be used if the working reagent has a reagent blank absorbance less than 1.0 at 340 nm.

Precautions

1. This reagent is for *in vitro* diagnostic use only.
2. Avoid ingestion of reagent as toxicity has not yet been determined.
3. Reagents contain sodium azide as preservative

SPECIMEN COLLECTION AND STORAGE

1. Serum is recommended.
2. Plasma containing anticoagulants should not be used.
3. All material coming in contact with the sample must be free of ammonia and heavy metals.
4. Urea in serum is reported stable for 72 hours refrigerated at 2-8°C. Unrefrigerated sera should be used within eight hours.
5. All blood samples should be considered potentially infectious.

INTERFERENCES

1. Urease action is inhibited by fluoride.
2. Samples with abnormal ammonia levels give falsely elevated BUN results.
3. Bilirubin: No interference up to 40 mg/dL.
4. Hemoglobin: No interference up to 500 mg/dL.
5. Lipemia: No interference in the presence of triglycerides up to 1500 mg/dL.
6. Ascorbic acid: No interference up to 25 mg/dL.

MATERIALS REQUIRED BUT NOT PROVIDED

1. Accurate pipetting devices.
2. Timer.
3. Test tubes
4. Spectrophotometer with a temperature controlled cuvette able to measure at 340nm

PROCEDURE

Wavelength	: 340 nm
Temperature	: 37°C
Optical path	: 1 cm
Assay Type	: Fixed time
Reaction Direction	: Decreasing

1. Bring to room temperature (15 -30 °C)
2. Set the photometer to 0 (zero) with distilled water.

	Standard	Sample
Reagent 1	800 μ L	800 μ L
Standard	10 μ L	--
Sample	--	10 μ L

1. Mix and wait one minute.

	Standard	Sample
Reagent 2	200 μ L	200 μ L

2. After 30 sec. read and record the absorbance(Abs.1)
3. Return tube to 37 °C and after 60 sec.read the final Abs.(Abs.2)
4. Determine the Δ Abs for sample and standard

CALIBRATION

Use a UREA standard or serum calibrator. The procedure should be calibrated according to the instrument manufacturer's calibration instructions. If control results are found to be out of range, the procedure should be recalibrated.

CALCULATION

(A₁-A₂) = Absorbance change

(A₁ - A₂) Sample x concentration = UREA (mg/dL)

(A₁ - A₂) standard of standard

CONVERSION FACTOR: mg/dL x 0.1665 = mmol/L.

QUALITY CONTROL

The validity of the reaction should be monitored by use of the control sera with known normal and abnormal UREA values.

EXPECTED VALUES

10 - 50 mg/dL

It is strongly recommended that each laboratory establish its own reference range.

PERFORMANCE

Sensitivity/ Limit of Detection(LOD):

The lower limit of detection is 1,17 mg/dL.

Linearity:

When the recommendation is designed and studied, the test is linear up to 250 mg/dL. Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

Precision :

	Within run (n=20)		Between run (n=20)	
Mean (mg/dL)	37,5	110,6	37,0	110,3
SD	0,64	1,34	0,56	1,66
CV(%)	1,71	1,21	1,51	1,51

Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x).The results obtained using 40 samples were the following:

Correlation Coefficient (r) =0,9985

Regression equation y=1,0023x +0,157

The results of the performance characteristics depend on the analyzer used.

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BİLİMSEL TIBBİ ÜRÜNLER PAZ.SAN.VE TİC.LTD.ŞTİ