
REF	CONT	REF	CONT	REF	CONT
URI-10100	2 x 50 mL	URI-10180M	4 x 45 mL	URI-10100A2	5 x 20 mL
URI-10300	6 x 50 mL	URI-10180P	4 x 45 mL		
URI-10600	6 x 100 mL	URI-10240M2	6 x 40 mL		
URI-10180A	4 x 45 mL	URI-10240M3	6 x 40 mL		

**KULLANIM AMACI**

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
INTENDED USE


For the *in vitro* quantitative determination of Uric Acid in serum.

METHODOLOGY

This method uses uricase, peroxidase and TOOS the chromogen to yield a colorimetric end product. The colorimetric end product produced in this reaction can be measured at 520 nm and is proportional to the uric acid concentration in the sample.

Principle

Uric Acid present in the sample originates a coloured complex according to the following reactions;

REAGENT COMPOSITION

Active Ingredients

	Concentration
TOOS	<1 mM
4-AAP	<4 mM
Uricase	150 U/L
Peroxidase	12000 U/L
Non-reactive stabilizers and fillers.	
pH 7.6±0.2	

Precautions:

- This reagent is for *in vitro* diagnostic use only.
- Reagent contains Sodium Azide as a preservative.
- Serum specimen should be considered infectious and handled appropriately.

REAGENT PREPARATION

Reagent is ready to use form.

REAGENT STORAGE

- Store the reagent at 2-8°C (refrigerated).
- The reagent is stable until the expiration date when stored at 2-8°C.

REAGENT DETERIORATION

Do not use the reagent if:

- The reagent is turbid.
- The reagent blank has an absorbance of 0,3 or greater at 520nm.
- The reagent does not meet stated performance parameters.

SPECIMEN COLLECTION AND STORAGE

- Test specimen should be serum and free from hemolysis.
- Bacterial contamination should be avoided to preserve the loss of uric acid.
- Uric acid in serum is stable for 3 days at 2-8°C, and up to 6 month when frozen.

INTERFERENCE

- Bilirubin: No interference up to 40 mg/dL.
- Lipemia: No interference in the presence of triglycerides up to 1500 mg/dL.
- Collection tubes containing formaldehyde as a preservative must be avoided.

ADDITIONAL EQUIPMENT REQUIRED BUT NOT PROVIDED

- A clinical chemistry analyzer capable maintaining constant temperature (37°C), and measuring absorbance at 520 nm.
- Deionized water and related equipment, e.g.: pipettes
- Analyzer specific consumables, e.g.: sample and read cups.
- Control and calibrator materials.

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Rev.Date / No: 30.08.2023 / 10

	Blank	Standard	Sample
Reagent	1000 µL	1000 µL	1000 µL
Distilled water	25 µL	--	--
Standard	--	25 µL	--
Sample	--	--	25µL

Mix and then incubate for 10 min at 37°C. Measure the absorbance of sample and standard against the reagent blank. The colour is stable at least 30 min.

CALCULATIONS

(Abs = Absorbance)

$\frac{\text{Abs (Sample)}}{\text{Abs (Standard)}} \times \text{Concentration of standard (mg/dL)} = \text{Uric Acid (mg/dL)}$

LIMITATIONS:

The reagent and sample volumes may be altered proportionally to accommodate various instrument requirements.

CALIBRATION

Use an aqueous Uric Acid standard or serum calibrator.

QUALITY CONTROL

The integrity of the reaction should be monitored by use of a two level control with known Uric Acid values. If control values are found outside the defined range, check the instrument, reagents and technique for problems.

EXPECTED VALUES

1,5 – 7,0 mg/dL

It is strongly recommended that each laboratory establish its own normal range.

PERFORMANCE

Sensitivity/ Limit of Detection(LOD):

The lower limit of detection is 0,063 mg/dL.

Linearity:

When the recommendation is designed and studied, the test is linear up to 25 mg/dL. Samples exceeding this value are diluted 1+1 with NaCl (9 g/L) and the result is multiplied by 2.

Precision :

	Within run (n=20)		Between run (n=20)	
Mean (mg/dL)	5,02	11,49	4,98	11,15
SD	0,05	0,10	0,03	0,10
CV(%)	1,04	0,84	0,69	0,91

Accuracy (Method Comparison) :

Results obtained BTPRODUCTS Reagents (y), did not show systematic differences when compared with other commercial reagents(x).The results obtained using 40 samples were the following:

Correlation Coefficient (r) =0,9973

Regression equation $y = 0,8774x + 0,3307$

The results of the performance characteristics depend on the analyzer used.

REFERENCES

- Davidson, L.and Henty,J.B.:Todd Sanford Clinical Diagnosis by laboratory Method,15 th ed.W.B.Saunders Company,Philadelphia,PA(1974).: Diagnostic
- Caraway, W.T., Clin. Chem. 4:239 (1963).
- Morin, L.G., Clin. Chem. 20:51 (1974).
- Fossati,P.,Principe L.,and Bertia,A.,Clin.Chem.26:227(1980)
- Duncan, p. Et al.,Clin.Chem.28:291(1982).
- Henry R.,J.Clin.Chem.:Principles and Techniques NY,Harper and Row,Second ed. (1974).
- Young, D.S., et al. Clin. Chem. 21:10-4320 (1975).
- Tietz,N.W.,Fundamentals of Clinical Chemistry,Philadelphia,W.B.Saunders729(1976)