KULLANIM AMACI

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INTENDED USE

For the quantitative determination of zinc in serum ,plasma or urine.
For in vitro diagnostic use only.

METHOD HISTORY

5-BR-PAPS use in routine zinc determinations because, this compound is water soluble, and has a higher molar absorptivity. Results by the proposed method developed with 5-Br-PAPS correlated well with those by atomic absorption spectrophotometry. This method is fast,sensitive and linear over the entire concentration range of

zinc in human serum, plasma or urine.

PRINCIPLE

Zinc produces with the specific complexant 5-Br-PAPS [(2-5-bromo-2-pyridylazo)-5-(N-propyl-N-sulfo-propylamino) phenol] a stable coloured complex, which colour intensity is proportional to the amount of zinc in the sample. The interference due do the oligoelements (iron and copper) present in the sample, are eliminated using peculiar reaction conditions and specific masking agents.

REAGENT COMPOSITION**Active Ingredients**

	Concentration
Bicarbonate buffer ph 9. 8	185 mM
5 Br-PAPS	0.05 mM
Sodium citrate	125 mM
Dimethylglyoxime	5 mM
Surfactants	% 1
Preservatives and stabilizers	

REAGENT PREPARATION

Reagents are ready to use .

STORAGE AND STABILITY

When stored at 15-25° C, the components of the kit will remain stable until the expiration date stated on the label.

PRECAUTIONS

1. The reagents contain Sodium azide.
2. Handle with care. The disposal of the residues has to be made according to legal local regulations.
3. To avoid ingestion and contact with eyes, skin and mucous membranes.

SPECIMEN COLLECTION AND STORAGE

1. Serum or nonhemolyzed plasma:Use only heparin salts as anticoagulant.
2. 24 hour urine : Diluted 1:10
3. Shake and bring the samples at room temperature (15 – 25°C) before using.
4. Turbid and icteric sera can interfere in the reaction raising values.

INTERFERENCES

EDTA anticoagulant masks zinc to 5-Br-PAPS chromogenic system. the test is not affected by the presence of conjugated and non-conjugated bilirubin up to 15 mg/dL, hemoglobin up to 0.5 g/dL and triglycerides up to1000 mg/dL, copper up to 1000µg/dL,iron up to 1000µg/dL and cobalt up to 500 µg/dL.

MATERIALS REQUIRED BUT NOT PROVIDED

1. Clinical chemistry analyzer capable maintaining constant temperature (37°C) and measuring absorbans at (560)nm
2. Deionize water and related equipment e.g. pipettes
3. Analyzer specific consumables, e.g: sample and read cups
4. Control and calibrator materials

PROCEDURE

Wavelength	: 560 nm
Working temperature	: 37°C
Optical path	: 1 cm
Assay type	: Endpoint
Direction	: Increasing

	Blank	Standard	Sample
Reagent	1000 µL	1000 µL	1000 µL
Standard	--	50 µL	--
Sample	--	--	50 µL

Mix well and let stand 10 min. at room temperature (20-25°C), or 5 min.37°C.And measure the absorbance of sample and standard against blank .

CALCULATION

(Abs=absorbance)

Sample Abs. x std.conc. = ZINC (µg/dL)

Standard Abs.

Conversion factor : (µg/dL) x 0.153 =µmol/L

In the urine : multiply result VOL.(dL) of 24 h-urine by 10(dilution factor) .

CALIBRATION

bt products ZINC standard must be used.

QUALITY CONTROL

The validity of the reaction should be monitored by use of the control sera with known normal and abnormal ZINC values.

EXPECTED VALUES

Men : 60 – 110 µg/dL

Women : 60 – 107 µg/dL

Urine : 150 – 1200 µg/24 h.

These values should only be used as a guideline.

Each laboratory should establish its own reference range.

PERFORMANCE

1. **Linearity limit:** The reaction is linear to 400 µg/dL.For higher value dilute sample 1:3 with deionized water , repeat test and multiply the result by 3 .
2. **Detection Limit** : 5 µg/dL . The lower detection limit represents the lowest measurable Zinc concentration that can be distinguished from zero.
3. **Precision:**

	Within run			Run to Run		
Mean(µg/dL)	77.55	202.65	300.65	77.55	202.65	300.65
SD	2.29	4.60	3.38	1.14	0.90	3.27
CV(%)	2.95	2.27	1.13	1.46	0.45	1.09

4. **Accuracy:** Comparison of Zinc colorimetric (y) with an other commercial reagent (x) assay gives the following result:

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Bilimsel Tıbbi Ürünler Paz.San.Tic.Ltd.Şti.