
REF
CWS-201000
CWS-20180A

Packing Size/ Ambalaj	CONT
1 x 1000 mL	1 x 1000 mL
4 x 45 mL	4 x 45 mL

**CUVETTES' WASHING SOLUTION**

The cuvettes' washing solution is ready to use, exactly as it is presented. It is suitable for washing cuvettes, flow cells and hydraulic circuits of photometer and autoanalyzers.

USE INSTRUCTIONS

CWS-201000 : For only use Automatic Analyzers and photometers

BIOTECNICA INSTRUMENTS

BT 1000 / 2000 / 3000 / 3500 / 35i

Proceed daily at the end of the working day with one of the following procedures.

1. Automatic procedure: at the end of the working day switch off the instrument with the "Shut Down" procedure and follow the instructions given by the program.

2. Upon request procedure:

BT 3000 TARGA :

Enter in the "System" menu, and then select the "Diagnostic" line; in the "diags" page select the "Function" menu and then the "Totally wash" command. Follow the instructions given by the program.

BT 35i / 3500 / BT 3000 PLUS / 2000 PLUS / 1000 :

Click on the "Extra Wash" button in the "Analyzer Utilities" window and follow instructions.

Bear in mind that the analyzer will use about 15 mL of solution for washing cycle.

ALL MANUAL PHOTOMETERS

Let aspirate about 10 mL of washing solution, then allow abundant rinsing with distilled water.

ATTENTION: The cuvettes' washing solution must be kept tightly closed and stored at 15-25°C.

H314 Causes severe skin burns and serious eye damage.
P303+P361+P353 IF ON SKIN:Remove/take off immediately all contaminated clothing.Rinse skin with water/shower
P304+P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
P305+P351+P338 IF IN EYES :Rinse continuously with water for several minutes .Remove contact lenses if present and easy to do.Continue rinsing

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Bilimsel Tıbbi Ürünler Paz.San.Tic.Ltd.Şti.