REF	Packing Size/ Ambalaj	CONT
EWS-10BUL	1 x 1000 mL	1 x 1000 mL
EWS-10180A	4 x 45 mL	4 x 45 mL

: 15-25 °C  **CUVETTE'S EXTRA WASHING SOLUTION**

The cuvettes' extra washing solution is ready to use, exactly as it is presented. It is suitable for washing cuvettes, flow cells and hydraulic circuits of photometer and autoanalyzers.

**REAGENT**

HCL and surfactant .

**USE INSTRUCTIONS**

**EWS-10** : For only use Automatic Analyzers and photometers  
**BIOTECNICA INSTRUMENTS**  
**BT 1000 / 2000 / 3000 / 3500 / 35i**

Proceed daily at the end of the working day with one of the following procedures.

**1. Automatic procedure:** at the end of the working day switch off the instrument with the "Shut Down" procedure and follow the instructions given by the program.

**2. Upon request procedure:**

**BT 3000 TARGA :**

Enter in the "System" menu, and then select the "Diagnostic" line; in the "diags" page select the "Function" menu and then the "Totally wash" command. Follow the instructions given by the program.

**BT 35i / 3500 / BT 3000 PLUS / 2000 PLUS / 1000 :**

Click on the "Extra Wash" button in the " Analyzer Utilities" window and follow instructions.

Bear in mind that the analyzer will use about 15 mL of solution for washing cycle.

**ALL MANUAL PHOTOMETERS**

Let aspirate about 10 mL of washing solution, then allow abundant rinsing with distilled water.

**ATTENTION: The cuvettes' extra washing solution must be kept tightly closed and stored at 15-25°C.**

H314 Causes severe skin burns and serious eye damage.  
H335 May cause respiratory irritation  
P303+P361+P353 IF ON SKIN:Remove/take off immediately all contaminated clothing.Rinse skin with water/shower  
P304+P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.  
P305+P351+P338 IF IN EYES :Rinse continuously with water for several minutes .Remove contact lenses if present and easy to do.Continue rinsing

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Bilimsel Tıbbi Ürünler Paz.San.Tic.Ltd.Şti.