REF

CA21-40

CONT

1*40 mL

: 2-8 °C  IVD

Please carefully read this package insert before use.

Only use for In Vitro Diagnostic.

INTENDED USE

For calibrating WBC and RBC and monitoring or evaluating the precision of test results on urine sediment analyzer.

PRINCIPLE

This product is colorless to pale yellow cell suspension. Base on flow cell digital imaging principle, calibrator is sampled into the test channel for picture collection and counting, so as to get the test report.

INGREDIENT

Veterinary erythrocytes, anticoagulant, preservative solution, cell fixative and etc.

STORAGE AND STABILITY

- Storage: Stored at 2-8 °C in a sealed container. Avoid freeze and strenuous vibration.
- Stability: Product could be kept for 3 months if stored sealed at 2-8 °C. Once opened, period of validity is 30 day.

INSTRUMENT

This product are used on BT URICELL 1280 and URIT UD 1320 Urine Sediment Analyzers.

PROCEDURE

- Remove a vial of calibrator from refrigerator and leave it warm to room temperature (15-30 °C) for 3 to 5 minutes.
- Gently invert the vial for 1 to 2 minutes at the first time. For subsequent use, invert 8 to 10 times again is enough.
- Check whether the calibrator has been mixed well. If not, repeat STEP 2 until completely mixed.
- Take out the calibrator for routine test or quality control according to the requirements in Urine Sediment Analyzer Operation Manual.
- Compare the test result with the reference value offered by manufacturer.
- Tightly cap the vial and store it in refrigerator after use.

EXPLANATION VALUES

Mean values and acceptable ranges are different from lot to lot. Please refer to the Target Value Sheet attached to the product. It is recommended that each laboratory should establish its own mean values and acceptable ranges. Use the values provided in Operation Manual as reference only.

LIMITATION

- This product is used for monitoring BT URICELL 1280 and BT URIT UD 1320 Urine Sediment Analyzers.
- Inability to obtain expected values indicates product deterioration. Deterioration of the product may be caused by overheating or freezing during shipping or storage.

PRODUCT DETERIORATION

- Do not follow the specified condition and procedures to store and use this product may cause product deterioration.
- Deteriorated product will cause results errors. Dispose of it as waste blood.

NOTICE

- Use this product under the professional technician's guidance.
- Read the instruction carefully before use.
- Store as required.
- Ensure that the reagents match with the analyzer model, or

improper reagents may cause results errors.

- Inactivation disposal has been done, but protection measures are also needed during use.

CAUTION

- Dispose of the product as waste blood if opened for 24 hours or over the shelf life.
- Disposal method of waste blood: Add 30mL sodium hypochlorite solution with 5% effective chlorine into 1L waste solution, mix and keep them for 30 minutes, then drain out and flush with plenty of water.
- Blood cell deposition will not affect the values if within shelf life.

PERFORMANCE**1. Homogeneity****1.1. In-vial homogeneity**

The In-vial homogeneity should be complied with table 1:

Table 1

Parameter	RBC	WBC
CV / %	≤15.0	≤15.0

1.2. Vial-to-vial homogeneity

The vial-to-vial homogeneity should be complied with table 2:

Table 2

Parameter	RBC	WBC
CV / %	≤15.0	≤15.0

2. Biological Safety

Test results of HBsAg ,HIV-1/HIV-2 antibody ,HCV antibody are all negative

Lot Number	1821426	
	Unit	Reference Value
RBC	Cell / uL	202 ± 61
WBC	Cell / uL	201 ± 60

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Bilimsel Tıbbi Ürünler Paz.San.Tic.Ltd.Şti.