**INTENDED USE**

SHEATH Reagent is used in urine sediment analyzers (Urit UD 1320 and others) to dilute the sample and coat cells to be able to count and classify cells.

**FUNCTIONS**

During testing, the analyzer uses the urine sample screened by Sheath Reagent, creating a stable and horizontal flow through the analyzer. Sample images are captured and divided by the optical imaging system. Finally, the analysis results are defined by artificial neural network recognition technology.

**REAGENT COMPOSITION**

Buffer and detergent

**REAGENT PREPARATION**

Ready to use reagent without preparation.

**STABILITY AND STORAGE**

- 1.Store in dark place ,keep away from direct sunlight ,temperature should be 2 - 35°C.
- 2.Expire date 24 month that sealed storage and avoid exposing to direct sunlight ; 60 days after unscrewing the detergent container.

**OPERATIONAL PROCEDURE**

Details please refer to the operation manual of the urine sediment analyzer.

**SPESIFICATIONS**

- 1.The pH value is 7,6±0,2 ( 25 °C ± 1 °C )
- 2.Blank value: when you do the blank test in the urine sediment analyzer ,the total number of particles value should be 0.

**RESTRICTION,PRECAUTION AND NOTICE**

- 1.Use fresh human urine sample.The polluted urine sample is unavailable
- 2.Use the Sheath under professional's supervision. It is for in vitro diagnosis only.
- 3.Do not drink. In case of accidental drinking, seek medical advice.
- 4.Avoid contact with skin and eyes. In case of contact, rinse immediately with plenty of water or see a doctor.
- 5.Dispose the waste, rest products and polluted wrappage in accordance with local regular's requirements

**REFERENCE RANGES**

	RBC	WBC	CAST
Male	0-10/µL	0-15/µL	0-1/µL
Female	0-15/µL	0-15/µL	0-1/µL

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Bilimsel Tıbbi Ürünler Paz.San.Tic.Ltd.Şti.